

Fuels, Engines, and Emissions

Characterization of Transient Process Chemistry via Advanced Diagnostics

Background

High-efficiency advanced diesel engine technology is a leading near-term option for reducing petroleum consumption and oil imports in the United States. However, new regulations by the U.S. Environmental Protection Agency (EPA), to be phased in over the 2007–2009 model years, will require a 90% reduction in emissions of nitrogen oxides (NO_x) and particulate matter. Without development of improved engine-control strategies, new emission control systems and more advanced fuels, it is unlikely that diesel engines will achieve the new emissions standards.

Likely diesel solutions will combine advanced fuels, combustion regimes and catalysts. Catalysts for NO_x , sulfur and particulate control, are typically applied to open- or closed-cell honeycomb type structures which contain small channels millimeters wide and several inches long. In operation, these engine-catalyst systems exhibit temporally and spatially varying chemistry distributions throughout the

system: e.g., engine-generated hydrogen (H_2) pulses to effect catalyst performance. Measurement strategies capable of resolving such dynamic species distributions, including the intra-catalyst-channel chemistry variations, are required to develop efficient diesel solutions to the emissions challenge.


The Technology

The spatially resolved capillary-inlet mass spectrometer (SpaciMS) is a transportable instrument for quantifying the spatial and temporal distributions of emissions. The instrument was developed specifically for diesel catalysis research, but is applicable to other areas of research as well.

The SpaciMS uses a minimally invasive capillary inlet system (200-micrometer outer diameter, 10-milliliter per minute sampling rate) to transport time-varying species pools to a small mass spectrometer for analysis. A



Less dependence on foreign oil, and eventual transition to an emissions-free, petroleum-free vehicle



The SpaciMS output shows the concentration of selected chemical species as a function of cycle time and location within the catalyst block.

vehicle systems

fuels & lubricants

emission control

multiport valve selects between up to twelve capillary sampling locations deployed throughout a test system (e.g., an engine-catalyst system). The minimally invasive nature of the sample capillaries allows for in situ analysis of full intra-catalyst spatiotemporal species distributions.

With respect to understanding detailed catalyst chemistry this is a quantum leap from the previously available technology, which provided only catalyst-in and catalyst-out measurements, at much lower speeds. The broad species applicability of the SpaciMS allows analysis of a range of species critical to diesel catalysis, including NO_x, H₂, oxygen (O₂), carbon dioxide (CO₂), hydrocarbon fragments, sulfur dioxide, and hydrogen sulfide. The readily transportable nature of the SpaciMS has allowed for field application at industrial research laboratories.

Commercialization

The SpaciMS has been applied to make the first known measurements of H₂ generated via advanced in-cylinder diesel combustion strategies. Hydrogen plays a key role in catalyst performance and was also tracked throughout a catalyst system.

The SpaciMS has been used to monitor H₂ generated via the water-gas-shift reaction; the reaction of carbon monoxide (CO) and water (H₂O) over platinum to generate H₂ and CO₂. The figure on the previous page shows the intra-catalyst dynamics of H₂, CO₂

and H₂O associated with this reaction throughout a 3-inch-long catalyst. The phase and amplitude of the various dynamic species generations and depletion help to elucidate the detailed process chemistry.

ORNL has provided nine companies with information required to build and/or operate this system.

Benefits

- Broad species applicability
- Minimally invasive
- Quantifies intra-catalyst-channel chemistry
- High temporal resolution
- Easily transportable
- Measures undiluted diesel engine exhaust

Where Can I Find More Information?

Dr. Ronald L. Graves
Oak Ridge National Laboratory

865-946-1226

gravesrl@ornl.gov

<http://www.feerc.ornl.gov>

DOE Technology Manager
Gurpreet Singh

Department of Energy

202-586-2333

gurpreet.singh@hq.doe.gov

A Strong Energy Portfolio for a Strong America

Energy efficiency and clean, renewable energy will mean a stronger economy, a cleaner environment, and greater energy independence for America. Working with a wide array of state, community, industry, and university partners, the U.S. Department of Energy's Office of Energy Efficiency and Renewable Energy invests in a diverse portfolio of energy technologies.



U.S. Department of Energy

Energy Efficiency and Renewable Energy

December 2003