

# **ISOLEX™\* 300**

**Magnetic Cell Separator**

Yielding

*the*

Purity...



# ... Of Immunomagnetic Selection.





## Specific Identification.

The process begins with the 9C5 antibody, a highly specific, high-affinity murine anti-human CD34 antibody, manufactured under GMP quality standards by Baxter. Within minutes, the 9C5 antibody coats and identifies the target CD34+ cells, leaving other subpopulations unmarked.

## Specific Capture.


The marked stem cells are captured by rosetting them with magnetic beads coated with a sheep anti-murine (SAM) antibody. The SAM antibody recognizes only the murine 9C5 antibody labeling the target cells, resulting in highly specific capture of CD34+ cells. In addition very efficient capture of CD34+ cells is achieved by using small 4.5  $\mu\text{m}$  beads. The large total surface area that results (2-3m<sup>2</sup> per run) means an extremely large number of binding sites for capture. The rosetted CD34+ cells are then physically separated from the unrosetted cell population by means of the precisely-engineered Isolex™ magnet.



## Specific Release.

After the non-target cells are drained away, a pharmaceutical grade enzyme, Chymo-Cell-T™ (Chymopapain), specifically releases the captured CD34+ cells from the beads, allowing the CD34+ cells to be recovered with a high degree of efficiency. Their surfaces are free of residual primary antibody, which might otherwise cause immunogenic reactions or sterically hinder further antibody-mediated processing.

PBSC Pre Selection: 1.2% CD34+




Selected Cells: 95.0% CD34+



## High Purification.

By specific identification, capture and release, the Isolex™ system ensures both high yield and purity. This high level of CD34+ cell purity is illustrated by the FACS analysis to the left.

Consistent with this high level of purity, colony culture assays also show 50-200 fold enrichment of CFU-GM, BFUe and CFU-GEMM after Isolex™ selection.



The Isolex™ 300 Magnetic Cell Separator combines several carefully designed elements such as a specially engineered magnet array and an automated mixing device. User control is exerted via keypad and digital display run by an integrated software program, without requiring an additional stand-alone computer. The entire system is compact enough to be operated on the laboratory benchtop.

The primary magnet is integrated into the mixing device, and moves into position at the appropriate time to accomplish the separation. Automated separation occurs inside the 300 mL cylindrical chamber, the geometry of which enables gentle mixing and efficient purification. The secondary magnet and its disposable chamber are designed to trap the few beads which may have escaped the primary magnet. Both chambers are part of an integrated biocompatible disposable set that provides a closed aseptic flowpath.

## ISOLEX™ 300

Magnetic  
Cell Separator

ISOLEX 300  
MAGNETIC CELL SEPARATOR  
Baxter

## A. Background of Safety

Before initiating clinical trials, the Isolex™ 300 system and its component parts passed extensive preclinical safety testing.

Preclinical Tests <sup>10</sup>	Results
Mouse Acute Systemic Infection (Beads)	No Adverse Effects
Rabbit Intracutaneous Infection (Beads)	No Tissue Reaction
Colony Formation After Exaggerated Exposure to Magnetic Field Magnetic Beads CD34 Antibody Chymopapain	No Inhibition of Trilineage Colony Formation
Effect of Chymopapain on CD34 Antigen	βC5 Epitope Re-expression > 90% After 24 Hours; βG12 Epitope Unaffected by Chymopapain


## Documented GMP Quality

In accordance with GMP standards, a large battery of tests support the quality of the manufacturing process and the final products:

<b>Dynabeads®</b>  50 Separate Tests	<b>Key Process and Release Tests</b> Sterility Endotoxin Functionality Identity Particles Stability
<b>Isolex™ Anti-CD34 Antibody</b>  120 Separate Tests	<b>Key Process and Release Tests</b> Sterility* Endotoxin Functionality Identity Stability Purity


\*NOTE: Virus removal and inactivation are achieved by hollow-fiber filtration/column chromatography and Solvent/Detergent treatment, respectively.

## Isolex Process



# For Effective:

## Stem Cell Purification



▶ Isolox™ 300 gives high purity and yield independantly of the stem cell source. High performance is maintained even at lower starting CD34 percentages. ◀

## Tumor Cell Depletion

STEM CELL SOURCE	ORIGIN	TARGET CELLS	LOG DEPLETION
G-mPBPC <sup>2</sup>	Donor	B Cells (CD19)	3.0
G-mPBPC <sup>3</sup>	Donor	Lymphoma cell lines	5.3
BM <sup>4</sup>	Patient	Solid Tumors	2.5


▶ In various studies, the Isolox™ Magnetic Cell Separator has provided substantial depletion of tumor cells and other non-target cells. ◀

## T-Cell Reduction

STEM CELL SOURCE	ORIGIN	TARGET CELLS	LOG DEPLETION
G-mPBPC <sup>2</sup>	Donor	CD5	3.0
G-mPBPC <sup>2</sup>	Donor	CD3	3.3
G-mPBPC <sup>7</sup>	Donor	CD3	4.0
BM <sup>8</sup>	Patient	CD3	2.5

▶ T-Cell numbers are significantly reduced after positive selection with Isolox™ system. ◀

## Hematopoietic Reconstitution



▶ Initial clinical studies show no significant differences in engraftment between Isolox™-selected and unselected marrow. ◀

## Address and Order Information

Code	Description
R9722	Isolex™* 300 Magnetic Cell Separator
R9720	Isolex™* 300 Disposable Set
R9500	Isolex™* 300 Stem Cell Reagent Kit

For inquiries or placing an order please contact your local Baxter Customer Service listed below:

### Austria

Baxter Ges. mbh  
Richard-Strauss-Str. 33  
1232 Wien  
Tel (+43) 1 - 61 04 80  
Fax (+43) 1 - 61 55 36 3

### Belgium

Baxter S.A.  
Rue Cplonel Bourg 105B  
1140 Brussels  
Tel (+32) 2 - 70 21 88 6  
Fax (+32) 2 - 73 63 68 4

### Denmark

Baxter A/S  
Gydevang 43  
3450 Allerod  
Tel (+45) 48 - 14 01 44  
Fax (+45) 48 - 14 02 44

### Finland

Baxter Oy  
P.O. Box 46  
01621 Vantaa  
Tel (+358) 0 - 50 73 54 73  
Fax (+358) 0 - 87 81 03 8

### France

Baxter S.A.  
Avenue Louis-Pasteur  
Boite Postal 56  
78311 Maurepas Cedex  
Tel (+33) 1 - 34 61 52 43  
Fax (+33) 1 - 34 61 50 25

### Germany

Baxter Deutschland GmbH  
Edisonstrasse 3-4  
85716 Unterschleißheim  
Tel (+49) 89 - 31 70 12 72  
Fax (+49) 89 - 31 70 15 56

### Ireland

Baxter Healthcare Ltd.  
Unit 7, Deansgrange Industrial  
Estate  
Blackrock  
Co. Dublin  
Tel (+353) 12 - 89 27 11  
Fax (+353) 12 - 89 27 04

### Italy

Baxter S.p.A.  
Viale Tiziano N. 25  
00196 Rome  
Tel (+39) 6 - 32 49 12 43  
Fax (+39) 6 - 32 21 64 3

### Netherlands

Baxter B.V.  
Kobaltweg 50  
3542 CE Utrecht  
Tel (+31) 30 - 48 89 11  
Fax (+31) 30 - 41 17 55

### Norway

Baxter AS  
Gjerdrumsvei 10 B  
0486 Oslo 4  
Tel (+47) 22 - 18 41 01  
Fax (+47) 22 - 18 41 09

### Spain

Baxter S.A.  
s/Solsones, No. 2 bajo, local 7  
Edificio Muntadas  
Parque de Negocios  
"MAS BLAU"  
El Prat de Llobregat  
08820 Barcelona  
Tel (+34) 3 - 47 87 16 2  
Fax (+34) 3 - 47 87 10 9

### Sweden

Baxter Medical AB  
Isafjordsgatan 30B  
16494 Kista  
Tel (+46) 8 - 632 64 00  
Fax (+46) 8 - 752 01 12

### Switzerland

Baxter A.G.  
Industriestr. 31  
8305 Dietlikon  
Tel (+41) 1 - 83 33 36 0  
Fax (+41) 1 - 83 30 42 0

### United Kingdom

Baxter Healthcare Ltd.  
Wallingford Road  
Compton, Newbury  
Berkshire, RG16 0QW  
Tel (+44) 635 - 20 63 24  
Fax (+44) 635 - 20 61 26

## References

- 1 Tseng-Law J et al, *Exp Hematol* 22:20 (1994) and Data on File, Baxter Healthcare Corp.
- 2 Lane TA et al, *Blood* 85(1):275 (1995)
- 3 Marolleau JP et al, *Blood* 84(10):370 (1994) and Data on File, Baxter Healthcare Corp.
- 4 Traycoff CM et al, *Exp Hematol* 22:554 (1994)
- 5 Paulus U et al, *Exp Hematol* 22:15 (1994)
- 6 C Civin, Personal Communication; (1994) and data on File, Baxter Healthcare Corp. (data obtained using MY-10 anti-CD34 antibody)
- 7 Dreger P et al, *Exp Hematol* 23:L147 (1995)
- 8 Broun ER, Personal Communication; (1994) and data on File, Baxter Healthcare Corp.
- 9 - Broun ER et al, *Exp Hematol* 22:362 (1994)
- 10 Data on File, Baxter Healthcare Corp.

Note: All data has been obtained using the Isolex™ 50 and Isolex™ 300 Magnetic Cell Separation systems

Isolex is a trademark of Baxter International Inc. and is registered in the United States Patent and Trademark Office. European registrations are pending.

ChymoCell-T™ is a registered trademark of The Boots Company PLC.

Dynabeads is a registered trademark of Dynal A.S., Oslo, Norway.