SEP 2 9 2004 KOH 1866 - AxSYM®Testosterone Page 1 of 4

510(k) Summary

Introduction

According to the requirements of 21 CFR 807.92, the following summary of safety and effectiveness provides details to understand the basis for a determination of substantial equivalence.

Submitter name, address, contact

Axis-Shield Diagnostics Ltd. The Technology Park Dundee DD2 1XA Scotland, UK Tel : +44 1382 422000 Fax: +44 1382 422088

Contact Person:	Susan Leonard
Email address:	susan_leonard@uk.axis-shield.com
Date Prepared:	July 7 th , 2004

Device Name

Proprietary Name: Abbott AxSYM[®] Testosterone Microparticle Enzyme Immunoassay (MEIA) test

Common name: Microparticle Enzyme Immunoassay (MEIA) for the determination of Testosterone.

Classification name: Testosterone test system

Predicate Device

Substantial equivalence to the Roche Elecsys[®] Testosterone assay (K964889) is claimed.

K041866 - AxSYM[®]Testosterone Page 2 of 4

510(k) Summary

Device Description

1

Competition format. Microparticle Enzyme Immunoassay (MEIA) for use on Abbott AxSYM[®] system.

Assay procedure:

- Incubate the sample with the antibody-coated microparticles.
- Add testosterone alkaline phosphatase conjugate and incubate.
- Transfer to matrix cell.
- Wash to remove unbound substances.
- Add substrate.
- Measure fluorescent product.

Intended Use

A Microparticle Enzyme Immunoassay (MEIA) for the quantitative determination of testosterone in human serum and plasma on the AxSYM System.

Comparison to Predicate Device

The AxSYM[®]Testosterone immunoassay is substantially equivalent to the Roche Elecsys[®] Testosterone assay cleared under K964889. Both products are intended for use in the

quantitative determination of testosterone.

The following information is presented in the Premarket Notification [510(k)] for AxSYM[®] Testosterone to support a substantially equivalent determination.

Similarities:

- Intended use
- Operating principles (competitive assay format using antigen-antibody interactions)
- Sample type (serum and plasma)
- Quantitative assay
- For use on automated instruments



510(k) Summary

Comparison to Predicate Device (continued)

Differences:

1

Parameter	Submission Device	Predicate Device
	AxSYM [®] Testosterone	Elecsys [®] Testosterone
Technology Format	Microparticle Enzyme Immunoassay (MEIA).	Electrochemiluminescence Immunoassay (ECLIA).
Capture Antibody	Anti-testosterone (mouse, monoclonal) antibody coated microparticles.	Biotinylated monoclonal anti-testosterone antibody (mouse). Becomes bound to streptavidin coated microparticles.
Competition Antigen	Testosterone:Alkaline Phosphatase conjugate.	Testosterone derivative labelled with a ruthenium complex.
Substrate	4-Methylumbelliferyl phosphate.	Not applicable.
Assay End-Point	Fluorescence	Electrochemiluminescence
Quantitation	Results are determined from a standard calibration curve (0, 0.2, 1.0, 2.5, 7.0, 15.0 ng/mL) generated and stored on the instrument.	Results are determined from a calibration curve generated on the instrument via a two-point calibration and a master curve provided via the reagent barcode.
Standardisation	Concentration is standardised gravimetrically and confirmed by GCMS (Gas Chromatography Mass Spectrometry).	The method has been standardised via ID-GCMS (Isotope Dilution Gas Chromatography Mass Spectrometry).
Reagent stability	Up to stated expiration date when stored at 2-8°C. 2 weeks on-board the AxSYM.	Up to stated expiration date when stored at 2-8°C. 4 weeks on-board the Elecsys 1010. 8 weeks on-board the E170/Elecsys 2010.

__ - AxSYM[®]Testosterone Page 4 of 4

510(k) Summary

Comparison to Predicate Device (continued)

Revised

>

Performance Characteristics:

Parameter	Submission Device	Predicate Device
	AxSYM [®] Testosterone	Elecsys [®] Testosterone
Precision	Within-run CV of 3.9% to	Within-run CV of 0.9% to
	8.3%.	4.6%
	Total CV of ≤14% for	Total CV of 1.6 to 7.4 %
	samples with testosterone	from 0.24 to 7.01 ng/mL.
	concentrations ≥0.85ng/ml	
	and $\leq 10\%$ for samples	
	≥2.10ng/ml.	
Functional Sensitivity	≤0.2ng/ml	Not stated
Analytical Sensitivity	<u>≤0.1ng/ml</u>	0.02ng/ml
Analytical Specificity	25 potential cross-reactants	12 potential cross-
• <u> </u>	tested. <1% cross reactivity	reactants tested at
	observed at the	40ng/mL. Cross-reactivity
	concentrations stated in the	ranged from 0-10.4%.
	pack leaflet (dependent on	
	cross-reactant).	
Interfering Substances	<15% interference at:	No interference at:
	Bilirubin – 20mg/ dL	Bilirubin – 25mg/dL
	Haemoglobin – 750mg/dL	Haemoglobin – 1.0g/dL
	Protein –10g/dL	Not listed
	Triglycerides – 500mg/dL	Not listed
	SHBG – 90nmol/L	Not listed
	Not listed	Lipemia – 1500mg/dL
	Not listed	Biotin – 30ng/mL

DEPARTMENT OF HEALTH & HUMAN SERVICES



Public Health Service

Food and Drug Administration 2098 Gaither Road Rockville MD 20850

SEP 2 9 2004

Ms. Susan Leonard Regulatory Affairs Manager Axis-Shield Diagnostics Ltd. The Technology Park Dundee, Scotland United Kingdom DD2 1XA

Re: k041866

Trade/Device Name: Abbott AxSYM® Testosterone Microparticle Enzyme Immunoassay (MEIA)

Regulation Number: 21 CFR 862.1150 Regulation Name: Calibrator Regulatory Class: Class II Product Code: JIS, CDZ Dated: July 7, 2004 Received: July 9, 2004

Dear Ms. Leonard:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html.

Sincerely yours,

Kan M. Corper MS, DV.M.

Jean M. Cooper, MS, D.V.M. Director Division of Chemistry and Toxicology Office of *In Vitro* Diagnostic Device Evaluation and Safety Center for Devices and Radiological Health

Enclosure

Revised

Indications for Use

510(k) Number (if known): K041866

Device Name: Abbott AxSYM[®] Testosterone Microparticle Enzyme Immunoassay (MEIA)

Indications For Use:

AxSYM[®] Testosterone is a Microparticle Enzyme Immunoassay (MEIA) for the quantitative determination of testosterone in human serum and plasma on the AxSYM System.

Testosterone monitoring is used clinically to diagnose and differentiate endocrine disorders. In males, these include hypogonadism, testicular failure, infertility, hypopituitarism and hyperprolactinemia. In females, polycystic ovary syndrome, adrenal hyperplasia, infertility, hirsutism, amenorthea, obesity and virilization can cause changes in serum testosterone levels.

The AxSYM[®] Testosterone assay is used as an aid in the investigation of infertility in males and of hirsutism and virilization in females.

The AxSYM[®] Testosterone Standard Calibrators are for the standard calibration of the AxSYM system when used for the quantitative determination of testosterone in human serum and plasma.

Prescription Use V (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Page 1 of 1

7

Office of in Vitro Diagnostic Device Evaluation and Safety

K041866 510(k)