

At the conclusion of this talk, you will be able to.....

- List the major changes in the 2004 NCCLS tables (M100-S14)
- Describe new test/report recommendations for Staphylococcus spp. including:
 - testing for inducible clindamycin resistance

 use of cefoxitin disk test to detect oxacillinresistant staphylococci



At the conclusion of this talk, you will be able to.....(con't)

- Discuss disk diffusion testing of Stenotrophomonas maltophilia and Burkholderia cepacia
- Describe the new reference guide for QC testing frequency when various test components are modified



NCCLS Standards - 2004

♦M100-S14 Tables (2004)

.....to be used with text documents explaining how to perform the tests....


M2-A8 Disk Diffusion (2003)

M7-A6 MIC (2003)



Reference Terminology

-when I refer to
- M100 -- this means the new tables (M100-S14)
- M2 -- this means the disk diffusion method (described in M2-A8)
- M7 -- this means the MIC method (described in M7-A6)



Updated information in M100-S14			
Vol. 24 No. 1	M100-S14		
Updated Information in This Edition			
document. There are several important changes to th Subcommittee on Antimicrobial Susceptibility Testing	This document includes all of the tables from the NCCLS Disk Diffusion (M2) susceptibility testing document. There are several important changes to the tables that have resulted from meetings of the Subcommittee on Antimicrobial Susceptibility Testing during 2000. Included below is a summary of the changes in this document, which superiode the tables published in 2020 and in earlier years.		
Summary of Major Changes in This Docum	Summary of Major Changes in This Document		
	The list includes the "major" changes in this document. Other minor or editorial changes have been made to the general formatting and to some of the table footnotes.		
Additions/Changes	Additions/Changes		
Introduction to Tables:			
Organism <i>Versivia pestis</i> - Added along with antimicrob (M7; Warning Table)	Organism <i>Yerzivia pestir</i> - Added along with antimicrobial agents that must not be reported as susceptible (M7; Warning Table)		
Suggested Grouping of Antimicrobial Agents:			
Antimicrobial agents that should be considered for testin added (M7; Table 1B)	g and reporting on potential agents of bioterrorism		
Enterobacteriaceae:			

Major Changes

♦ Enterobacteriaceae

- More on Salmonella spp. and using nalidixic acid to screen for fluoroquinolone resistance
- Pseudomonas aeruginosa and other non-Enterobacteriaceae
 - Move levofloxacin from Test / Report Group "U" to Group "B"
 - Disk diffusion breakpoints for Stenotrophomonas maltophilia and Burkholderia cepacia

Major Changes (con't)

◆ Staphylococcus spp.

- Inducible clindamycin resistance testing/reporting
 Cefoxitin disk test for mecA
- ♦ Coagulase-negative staphylococci
 - More on mecA and oxacillin MIC results
 - More on reporting β-lactams on oxacillin susceptible isolates

♦ Enterococcus faecalis

Predicting imipenem susceptibility from ampicillin results

Major Changes (con't)

Quality Control

- Reference Guide for QC testing frequency
- QC ranges for *E. coli* ATCC 35218 and β-lactam / β-lactamase inhibitor combination agents when using Haemophilus Test Medium (HTM)

- Oritavancin QC ranges

- Staphylococcus aureus ATCC 29213
- Enterococcus faecalis ATCC 29212
- Streptococcus pneumoniae ATCC 49619

Major Changes (con't)


New breakpoints

- Gemifloxacin
 - Enterobacteriaceae (FDA-approved for Klebsiella pneumoniae)
 - · Haemophilus spp.
 - Streptococcus pneumoniae
- Telithromycin
 - Staphylococcus spp.
 - Haemophilus spp.
 - Streptococcus pneumoniae

Major Changes (con't)

+Additions to Table 1A


- Haemophilus spp.
 - Gemifloxacin (Group C)
- Streptococcus pneumoniae
 Gemifloxacin (Group (B)



Major Changes (con't)

Potential agents of bioterrorism

- Separate tables for this group of organisms
- Bacillus anthracis, Yersinia pestis, Burkholderia mallei, Burkholderia pseudomallei



Salmonella and Fluoroquinolones (FQ)

• "FQ-susceptible strains of Salmonella that test resistant to nalidixic acid may be associated with clinical failure or delayed response in FQ-treated patients with extraintestinal salmonellosis. Extraintestinal isolates of Salmonella should also be tested for resistance to nalidixic acid. For isolates that test susceptible to FQs and resistant to nalidixic acid, the physician should be informed that the isolate may not be eradicated by FQ treatment. A consultation with an infectious disease practitioner is recommended."

M100-S14 (M2, M7); Table 2A

Salmonella spp. (blood)		
	<u>MIC (μg/ml)</u>	
ampicillin	>32 R	
ciprofloxacin	≤0.25 S	
ceftriaxone	≤0.5 S	
trimeth-sulfa	>4/78 R	
Test nalidixic acid on extraintes ciprofloxacin MICs of 0.12–1.0 μg/n of 2.0 μg/ml is Intermediate and ≥4.	nl; a ciprofloxacin MIC	

	Salmo	onella and	Ciproflo	xacin
	CIP MIC (µg/ml)	NCCLS Interpretation	Likely Mutation	Nalidixic acid
	≤ 0.06	S	none	S
	0.12- 1	S	one	R*
	≥ 4	R	two	R
:		s with extraintestina Q therapy; use nalic step mutants.		
	Ū	-	Threifall et al. 2001. I Butt et al. 2003. EID 9	



Salmonella	spp. (blood)
------------	--------	--------

	<u>MIC (µg/ml)</u>
ampicillin	>32 R
ciprofloxacin	≤0.25 S*
ceftriaxone	≤0.5 S
trimeth-sulfa	>4/78 R

*....if nalidixic acid is resistant, add comment such as..... "This isolate demonstrates reduced susceptibility to fluoroquinolones. For some patients with extraintestinal Salmonella infections with such isolates, the isolates may not be eradicated by fluoroquinolone treatment. ID consult suggested."

Salmonella and Ciprofloxacin

• "Salmonella spp. isolated from sterile sites or from patients that have failed FQ therapy should be tested for the MIC of ciprofloxacin or susceptibility to nalidixic acid. Those isolates for which the ciprofloxacin MICs are ≥0.125 µg/ml or resistant to nalidixic acid should be considered to have reduced susceptibility to FQs and physicians should be warned that clinical failure or delayed response may be associated with FQ treatment of infections caused by these isolates".

Poutanen and Low. 2003. CMN 25:97

Salmonella and Reduced Ciprofloxacin Susceptibility (MIC ≥0.12 μg/ml)

 Isolates uncommon in USA (www.cdc.gov/narms/annual/2001)

- UK study (1999) 23% in S. typhi
 Mostly travelers from India and Pakistan Threlfall et al. 2001. EID. 7:448.
- Nalidixic acid screen study n=1010 Salmonella; 50 isolates w/ reduced ciprofloxacin susceptibility:
 - Sensitivity 100%
 - Specificity 87%
 - Hakanen et. al. 1999. JCM. 37:3572.

New Disk Diffusion Breakpoints

- Stenotrophomonas maltophilia
 - levofloxacin
 - minocycline
 - trimethoprim-sulfamethoxazole
- Burkholderia cepacia
 - ceftazidime
 - meropenem
 - minocycline
 - 35°C; ambient air; 20-24 h incubation

Excerpt from Table 2B (M2).... "Zone Diameter Interpretive Standards and Equivalent MIC Breakpoints for *P. aeruginosa, Acinetobacter* spp., *S. maltophilia,* and *B. cepacia*"

Agent	R	1	S	Comments
Ceftazidime	≤ 14	15-17	≥18	
	≤ 17	18-20	≥ 21	For B. cepacia
Minocycline	≤ 14	15-18	≥19	May be reported for <i>S.</i> maltophilia and <i>B.</i> cepacia also
Trimeth- sulfa	≤ 10	11-15	≥16	May be reported for S. maltophilia also




Relocation of Levofloxacin in Table 1 for *Pseudomonas aeruginosa* and Other Non-Enterobacteriaceae

...new M100-S14

Group B Primary Test	ciprofloxacin levofloxacin
Report Selectively	
ald M100 C12	
010 101100-513	
old M100-S13 Group U	levofloxacin or
	lomefloxacin or
Group U	

M100-S14 (M2, M7); Table 1 & 2B



Major Changes M100-S14 Staphylococcus

Staphylococcus spp.

 "Macrolide resistant isolates of S. aureus and coagulase-negative Staphylococcus spp. may have constitutive or inducible resistance to clindamycin [methylation of the 23S rRNA encoded by the erm gene also referred to as MLS_B (macrolide, lincosamide, and type B streptogramin) resistance] or may be resistant only to macrolides (efflux-mechanism encoded by the msrA gene)."

M100-S14 (M2, M7); Table 2C

Staphylococcus spp. Erythromycin / Clindamycin


Mechanism	Determinant	Erythro	Clinda
Efflux	msrA	R	S
Ribosome alteration	erm	R	S*
Ribosome alteration	erm	R	R constitutive

msrA = macrolide streptogramin resistance

erm = erythromycin ribosome methylase


*requires induction to show resistance








Optional reporting strateg	y
Staphylococcus	s aureus
erythromycin	R
oxacillin	R
penicillin	R
vancomycin	S
"Contact laboratory if results nee	



"D Test" - positive reaction


Inducible clindamycin resistance (*erm-*mediated)

...another example



"D Test" – negative reaction

NO induction (*msr*A-mediated erythromycin resistance)




"D Test" - positive reaction

Inducible clindamycin resistance (*erm-*mediated)


Routine purity plate:

 Streak 1/3 of plate for confluent growth

• Place 2 µg clindamycin disk 15 mm from edge of 15 µg erythromycin disk



"D Test" positive and optic	onal comment
Staphylococcu	is aureus
clindamycin	R
erythromycin	R
oxacillin	R
penicillin	R
vancomycin	S
"This <i>S. aureus</i> is presumed to detection of inducible clinda Clindamycin may still be effectiv	mycin resistance.



"D Test" negative and optional comment			
Staphylococcus aureus			
clindamycin	S		
erythromycin	R		
oxacillin	R		
penicillin	R		
vancomycin	S		
"This S. aureus DOES NOT demonstrate			
inducible clindamycin resistance in vitro."			

Inducible Clindamycin Resistance - Incidence

Varies considerably geographically

Community-associated MRSA

- Frequently erythromycin-R clindamycin-S
- Often msrA-mediated mechanism (NOT inducible)

+ USA report 2002

- 617 S. aureus erythromycin-R clindamycin-S
- 50% NOT inducible resistance
- Fiebelkorn et al. 2003. JCM. 41:4740.


Staphylococcus spp.

 "The results of disk diffusion tests using a 30 μg cefoxitin disk and alternate breakpoints (see box at end of this table) can be used to predict mecA –mediated resistance in staphylococci."

M100-S14 (M2, M7); Table 2C

Disk Diffusion Screen for *mec*Amediated Resistance in Staphylococci

- Perform standard disk diffusion test with cefoxitin (30 μg) disk
- Incubate 24 h; however, results may be reported after 18 h, if resistant
- Report results for OXACILLIN, not cefoxitin





Disk Diffusion Screen for <i>mecA</i> - mediated Resistance in Staphylococci (con't)				
Cefoxitin zone (mm)				
S. aureus	≤ 19*	≥ 20**		
CoNS	≤ 24 *	≥ 25**		
 * Report as oxacillin resistant ** Report as oxacillin susceptible CoNS, coagulase-negative staphylococci 				
		M100-S14 (M2, M7); Table 2C		

Staphylo	ococcus	s - Oxa	acillin
MIC (µg/ml):	Suga	l m 4	Dee
-	<u>Susc</u>	<u>Int</u>	<u>Res</u>
S. aureus	≤ 2	-	≥4
CoNS	≤ 0.25	-	≥ 0.5
DD (mm):			
	Res	Int	<u>Susc</u>
S. aureus	≤ 10	11-12	≥13
CoNS	< 17	-	> 18



Oxacillin Breakpoints Coagulase-Negative Staphylococci


- May overcall resistance for species other than S. epidermidis (e.g. S. lugdunensis, S. saprophyticus)
- For serious infections with CoNS other than S. epidermidis, testing for mecA or PBP 2a may be appropriate for strains having oxacillin MICs of 0.5 – 2 µg/ml or oxacillin zones ≤17 mm
- If mecA or PBP 2a negative, report as oxacillin susceptible

M100-S14 (M2, M7); Table 2C


Oxacillin Breakpoints Coagulase-Negative Staphylococci (con't)

- For oxacillin-resistant strains (including PBP 2a or mecA positive strains), report all β-lactams resistant
- For oxacillin-susceptible strains, report any β-lactams tested according to results generated

M100-S14 (M2, M7); Table 2C









Major Changes M100-S14 Enterococcus faecalis



Enterococcus faecalis


 "Ampicillin susceptibility can be used to predict imipenem susceptibility provided the species is confirmed to be *E. faecalis*".


M100-S14 (M2, M7); Table 2D

If MD requests imipenem results on E. faecalis				
<i>E. faecalis</i> (blood)				
ampicillin	S			
vancomycin	S			
gent synergy	R			
strep synergy	S			
"Ampicillin-susceptible <i>E. faecalis</i> are imipenem susceptible"				









Excerpt from: "Reference Guide to QC Testing Frequency"				
(for ATCC QC strains after 20-30 consecutive days of satisfactory daily testing)				
	No. of days of consecutive QC testing required]
MIC test modification	1	5	20 or 30	Comments
Use new shipment or lot number	x			
Use new manufacturer of MIC test			x	In addition, perform in- house validation studies
Software update that affects AST results		х		Monitor all drugs not just those implicated in software modification
M100-S14 (M2, M7); Table 3B				



Excerpt from: "Reference Guide to QC Testing Frequency" (con't)

- Note 1: Addition of any NEW antimicrobial agent requires 20 or 30 consecutive days of satisfactory testing, prior to use of this guide.
- Note 2: QC can be performed prior to or concurrent with testing patient isolates. Patient results can be reported for that day if QC results are within the acceptable limits.

M100-S14 (M2, M7); Table 3B

Excerpt from: "Reference Guide to QC Testing Frequency" (con't)

- Note 3: Manufacturers of commercial or inhouse prepared tests should follow their own internal procedures and applicable regulations.
- Note 4: For troubleshooting out-of-range results, refer to M2-A8 or M7-A6, QC section.
- Note 5: Broth, saline and/or water used to prepare an inoculum does not require routine QC.

M100-S14 (M2, M7); Table 3B

More Examples..... Application of Table 3C


(consecutive days of daily QC required)


Inoculum preparation, convert from:

- Visual to Prompt 20 or 30 days
- Visual to photometer 5 days
- Prompt to photometer 5 days
- Instrument / software
 - Repair of instrument 1 day (or more)
 - Replace hardware (e.g. reader-incubator) 20 or 30 days

Haemophilus Test Medium (HTM)

• Addition of disk diffusion and MIC QC ranges for *E. coli* ATCC 35218 (ßlactamase producing strain) and ßlactam / ß-lactamase inhibitor drugs





Agent	Drug class	Route of administration	FDA approved
Daptomycin (Cubicin)	lipopeptide	IV	yes
Gemifloxacin (Factive)	fluoroquinolone	PO	yes
Oritavancin	glycopeptide	IV	no
Telithromycin (Ketek)	ketolide	PO	no

Daptomycin

- In vitro activity against gram-positive bacteria including MRSA and VRE
- Mode of action
 - Bactericidal
 - Requires physiologic calcium

Daptomycin (con't)

- Susceptibility testing media requirements
 - Mueller-Hinton broth 50 mg/L calcium chloride
 - Mueller-Hinton agar 28 mg/L calcium chloride
- Currently, no NCCLS breakpoints; FDA breakpoints available in product literature
- NCCLS QC ranges available
- No resistance reported to date among *S. aureus* or Group A or B Streptococcus

Gemifloxacin


- Active against respiratory pathogens including Streptococcus pneumoniae, Haemophilus influenzae, Moraxella catarrhalis, Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila
- Inhibits DNA synthesis through inhibition of both DNA gyrase and topoisomerase IV (dual target)


Oritavancin

- Bactericidal in vitro
- Active against most gram-positive pathogens including VRE

Telithromycin


- Active against respiratory pathogens (S. pneumoniae, H. influenzae, M. catarrhalis, Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila)
- Active against other gram-positive bacteria that have inducible MLS_B mechanism of resistance and does not induce resistance
- Active against S. pneumoniae resistant to erythromycin and clarithromycin regardless of resistance mechanism







- ♦ Bacteria included in M100-S14
 - Bacillus anthracis
 - Burkholderia mallei
 - Burkholderia pseudomallei
 - Yersinia pestis
- New Tables
 - 1B (Antimicrobial agents to test/report)
 - 2K (MIC interpretive standards)



Recap of ... "Summary of Comments and Subcommittee Responses"

- ESBLs testing
 - Report cefepime as resistant for ESBL producers
 - There is limited data on use of cephamycins (e.g. cefoxitin, cefotetan) for treating infections caused by ESBL producers
 - If ESBL confirmatory test negative, report results as tested (do not override to resistant)
 - Currently, only *E. coli* and *Klebsiella* spp. are addressed in NCCLS ESBL testing rules

Recap of.. "Summary of Comments and Subcommittee Responses" (con't)

Other GNR

- There are no specific NCCLS recommendations for testing for inducible β-lactamases. To help detect resistance to 3rd generation cephalosporins resulting from selection of derepressed mutants, repeat testing after 3-4 days is suggested for *Enterobacter, Citrobacter, Serratia.*
- mecA and coagulase-negative staphylococci
 - Discussed above

Recap of.. "Summary of Comments and Subcommittee Responses" (con't)

Incubation temperature

– range 33-35°C

- oxacillin Staphylococcus spp., 33-35°C (not >35°C)
- ♦ Haemophilus spp.

– β -lactamase testing only would not detect BLNAR strains

- MIC testing –frequency of performing colony count to QC inoculum

 Perform at least quarterly
- QC of commercial McFarland standards
 Follow manufacturers recommendations

Some Issues Under Discussion by NCCLS

- Staphylococcus spp. re-evaluate moxifloxacin, gatifloxacin, levofloxacin, ciprofloxacin breakpoints
- Acinetobacter examine correlation of disk and MIC results for β-lactams and tetracyclines
- Development of new Guideline for testing bacteria not currently addressed in NCCLS AST standards (e.g. Corynebacterium, HACEK, etc.)

Enterobacteriaceae β-Lactam Breakpoints and ESBL Issues

- Re-evaluation of β-lactam breakpoints for Enterobacteriaceae
 - Example: cefotaxime
 - Current Susceptible at ≤ 8 µg/ml
 - Proposed Susceptible at ≤ 1 or $\leq 2 \mu g/ml$
 - Substantial data needed
 - Goal is to more accurately detect all β -lactamase and other β -lactam resistance mechanisms with revised breakpoints
- Changing breakpoints commercial systems project it will take 3 years ...much \$\$\$\$\$!

