

---

# Guidance for Industry

## Current Good Manufacturing Practice for Medical Gases

### *DRAFT GUIDANCE*

This guidance document is being distributed for comment purposes only.

Comments and suggestions regarding this draft document must be submitted within 120 days of publication in the *Federal Register* of the notice announcing the availability of the draft guidance. Submit comments to Docket Management Branch (HFA-305), Food and Drug Administration, 5630 Fishers Lane, rm.1061, Rockville, MD 20852. All comments must be identified with the docket number listed in the notice of availability that publishes in the *Federal Register*.

For questions regarding this draft document contact Duane Sylvia (301) 594-0095, [sylvriad@cder.fda.gov](mailto:sylvriad@cder.fda.gov).

**U.S. Department of Health and Human Services  
Food and Drug Administration  
Center for Drug Evaluation and Research (CDER)  
May 2003**

**Compliance**

# Guidance for Industry

## Current Good Manufacturing Practice for Medical Gases

*Additional copies are available from:*

*Office of Training and Communications  
Division of Communications Management  
Drug Information Branch, HFD-210  
Center for Drug Evaluation and Research  
Food and Drug Administration  
5600 Fishers Lane, Rockville, MD 20857  
(Phone 301-827-4573)  
Internet: <http://www.fda.gov/cder/guidance/index.htm>*

**U.S. Department of Health and Human Services  
Food and Drug Administration  
Center for Drug Evaluation and Research (CDER)  
May 2003**

**Compliance**

*Contains Nonbinding Recommendations*

*Draft — Not for Implementation*

**TABLE OF CONTENTS**

<b>I. INTRODUCTION.....</b>	<b>1</b>
<b>II. STATUTORY AND REGULATORY REQUIREMENTS.....</b>	<b>2</b>
<b>III. ORGANIZATION AND PERSONNEL .....</b>	<b>3</b>
<b>A. Responsibilities of the Quality Control Unit.....</b>	<b>3</b>
<b>B. Personnel Qualifications.....</b>	<b>3</b>
<b>C. Consultants.....</b>	<b>5</b>
<b>IV. BUILDINGS AND FACILITIES .....</b>	<b>5</b>
<b>A. Design and Construction .....</b>	<b>5</b>
<b>B. Security .....</b>	<b>6</b>
<b>V. EQUIPMENT .....</b>	<b>6</b>
<b>A. Equipment Cleaning and Maintenance.....</b>	<b>6</b>
<b>B. Equipment Calibration.....</b>	<b>7</b>
<b>C. Computerized Systems .....</b>	<b>7</b>
<b>VI. COMPONENTS, CONTAINERS, AND CLOSURES .....</b>	<b>8</b>
<b>A. General Recommendations .....</b>	<b>8</b>
1. <i>Prefill Inspections for Cylinders.....</i>	<i>9</i>
2. <i>Prefill Inspections for Cryogenic Home Containers.....</i>	<i>11</i>
3. <i>Dedication of Large Cryogenic Containers to Medical Use Only.....</i>	<i>11</i>
4. <i>Prefill Inspections of Large Cryogenic Containers.....</i>	<i>11</i>
5. <i>Prefill Inspections for Permanently Mounted Cryogenic Containers.....</i>	<i>12</i>
<b>B. Retesting of Containers .....</b>	<b>12</b>
<b>VII. PRODUCTION AND PROCESS CONTROLS.....</b>	<b>12</b>
<b>A. Written Procedures.....</b>	<b>12</b>
<b>B. Charge-in of Components .....</b>	<b>13</b>
1. <i>Temperature/Pressure Readings (Boyle's Law).....</i>	<i>13</i>
2. <i>Valve Assembly Leak Testing.....</i>	<i>14</i>
3. <i>Heat of Compression.....</i>	<i>14</i>
<b>C. Calculation of Yield .....</b>	<b>15</b>
<b>VIII. PACKAGING AND LABELING CONTROLS .....</b>	<b>15</b>
<b>A. Materials Examination and Usage .....</b>	<b>15</b>
<b>B. Labeling Control.....</b>	<b>16</b>
<b>C. Packaging and Labeling Operations .....</b>	<b>16</b>
<b>D. Drug Product Inspection .....</b>	<b>18</b>
<b>E. Expiration Dating .....</b>	<b>18</b>

***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

<b>IX. HOLDING AND DISTRIBUTION .....</b>	<b>19</b>
<b>A. Warehousing Procedures .....</b>	<b>19</b>
<b>B. Distribution Procedures and Recalls.....</b>	<b>19</b>
<b>X. LABORATORY CONTROLS.....</b>	<b>20</b>
<b>A. General Controls.....</b>	<b>20</b>
1. <i>Sampling Plan</i> .....	21
2. <i>USP Oxygen Monograph</i> .....	21
3. <i>Calibration of Instruments</i> .....	22
<b>B. Testing and Release for Distribution.....</b>	<b>22</b>
1. <i>Liquid-to-Liquid Filling; Oxygen Only</i> .....	23
2. <i>Liquid-To-Gas; Filling Large Cryogenic Containers</i> .....	24
3. <i>Liquefied (gas on top of liquid) Compressed Gas</i> .....	24
4. <i>Gas Mixtures</i> .....	25
5. <i>Liquid Nitrogen</i> .....	25
<b>C. Alternate Testing Methods.....</b>	<b>26</b>
<b>D. Stability Testing .....</b>	<b>26</b>
<b>E. Reserve Samples.....</b>	<b>27</b>
<b>XI. RECORDS AND REPORTS.....</b>	<b>27</b>
<b>A. Record Retention .....</b>	<b>27</b>
<b>B. Equipment Cleaning and Use Log.....</b>	<b>28</b>
<b>C. Component, Drug Product Container, Closure, and Labeling Records .....</b>	<b>28</b>
<b>D. Master Production and Control Records.....</b>	<b>28</b>
<b>E. Batch Production and Control Records.....</b>	<b>29</b>
<b>F. Production Record Review.....</b>	<b>30</b>
<b>G. Laboratory Records.....</b>	<b>30</b>
<b>H. Liquid Supply (Certificate of Analysis (COA)).....</b>	<b>32</b>
<b>I. Distribution Records.....</b>	<b>33</b>
<b>J. Complaint Files .....</b>	<b>33</b>
<b>K. Reporting Deaths and Injuries .....</b>	<b>34</b>
<b>XII. RETURNED AND SALVAGED DRUG PRODUCTS.....</b>	<b>35</b>
<b>A. Returned Drug Products .....</b>	<b>35</b>
<b>B. Drug Product Salvaging .....</b>	<b>35</b>
<b>XIII. AIR SEPARATION PLANTS OR UNITS (ASU) .....</b>	<b>35</b>
<b>XIV. STORAGE TANK INSTALLATIONS AT HEALTH CARE FACILITIES.....</b>	<b>35</b>
<b>XV. MEDICAL GAS MIX-UPS.....</b>	<b>36</b>

*Contains Nonbinding Recommendations*

*Draft — Not for Implementation*

<b>XVI. CARBON DIOXIDE AND HELIUM MANUFACTURERS AND WHOLESALE DISTRIBUTORS .....</b>	<b>37</b>
<b>XVII. EMERGENCY MEDICAL SERVICE (EMS).....</b>	<b>38</b>
<b>XVIII. GAS-TO-GAS ADAPTERS.....</b>	<b>39</b>
<b>XIX. ALTERNATIVE APPROACHES.....</b>	<b>39</b>
<b>ATTACHMENT: MEDICAL GAS MIX-UPS .....</b>	<b>40</b>
<b>GLOSSARY.....</b>	<b>43</b>

*Contains Nonbinding Recommendations*

*Draft — Not for Implementation*

# Guidance for Industry<sup>1</sup>

## Current Good Manufacturing Practice for Medical Gases

This draft guidance, when finalized, will represent the Food and Drug Administration's (FDA's) current thinking on this topic. It does not create or confer any rights for or on any person and does not operate to bind FDA or the public. You can use an alternative approach if the approach satisfies the requirements of the applicable statutes and regulations. If you want to discuss an alternative approach, contact the FDA staff responsible for implementing this guidance. If you cannot identify the appropriate FDA staff, call the appropriate number listed on the title page of this guidance.

### I. INTRODUCTION

This guidance is intended to provide recommendations on how to comply with the current good manufacturing practice (CGMP) regulations as they apply to manufacturing, filling, transfilling, cascading, transferring, and distributing compressed and cryogenic medical gases. The recommendations should help manufacturers, fillers, and distributors<sup>2</sup> comply with CGMP requirements to ensure the identity, strength, quality, and purity of medical gases. This guidance also provides recommendations to medical gas manufacturers on how to comply with certain aspects of the PDMA final rule (i.e., 21 CFR part 205). This guidance is not intended to be an all-inclusive listing of all relevant CGMP; instead, it covers certain sections of the CGMP regulations followed by a discussion of recommendations that the Agency considers acceptable means of meeting the requirements.

Three previous documents were published on current good manufacturing practice for medical gases. FDA's first guideline on compressed medical gases was issued in June of 1981 and revised in 1983. In February of 1989, FDA issued another revision of the guideline to address the evolving home care area, including the delivery of oxygen to patients at home. This guidance builds on the previous guidelines. It provides details on the filling of high-pressure cylinders and cryogenic containers and includes new information on CGMP policy for large cryogenic containers, as well as discussion of CGMP relating to storage tank installation, carbon dioxide and helium manufacturing, and emergency medical services. Once finalized, this version of the guidance will supersede those earlier guidelines.

---

<sup>1</sup> This guidance has been prepared by the Division of Manufacturing and Product Quality in the Office of Compliance of the Center for Drug Evaluation and Research (CDER) at the Food and Drug Administration.

<sup>2</sup> For the purposes of this document, the term *manufacturer* includes fillers, transfillers, cascaders, distributors, and transferers of medical gases.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

39 FDA's guidance documents, including this guidance, do not establish legally enforceable  
40 responsibilities. Instead, guidances describe the Agency's current thinking on a topic and should  
41 be viewed only as recommendations, unless specific regulatory or statutory requirements are  
42 cited. The use of the word *should* in Agency guidances means that something is suggested or  
43 recommended, but not required.  
44

45

## **II. STATUTORY AND REGULATORY REQUIREMENTS**

47

48 Medical gases (e.g., oxygen, carbon dioxide, helium, nitrogen, nitrous oxide, medical air, and  
49 combinations of these) are drugs within the meaning of section 201(g)(1) of the Federal Food,  
50 Drug, and Cosmetic Act (the Act) (21 U.S.C. 321(g)(1)) and pursuant to section 503(b)(1)(A) of  
51 the Act (21 U.S.C. 353(b)(1)(A)) are required to be dispensed by prescription.

52

53 Medical gases are considered adulterated under section 501(a)(2)(B) of the Act (21 U.S.C.  
54 351(a)(2)(B)) if the methods used in, or the facilities or controls used during their manufacture,  
55 processing, packing, or holding do not conform to, or are not operated or administered in  
56 conformity with CGMP. The CGMP regulations are intended to ensure that a drug meets the  
57 safety requirements of the Act and has the identity and strength and meets the quality and purity  
58 characteristics that it purports or is represented to possess. Medical gases are finished drug  
59 products and are subject to the CGMP regulations at 21 CFR parts 210 and 211. Manufacturers  
60 of medical gases must follow the requirements in the CGMP regulations to comply with section  
61 501(a)(2)(B). For example, each time a medical gas is filled into another container, finished  
62 product testing must be performed in accordance with § 211.165(a).

63

64 Medical gases that are not produced and handled in accordance with CGMP regulations can  
65 cause serious injury or death to the patients who use them. A number of injuries and deaths have  
66 resulted from mix-ups of medical gases associated with CGMP violations including:

67

- 68 • Mislabeling (in some cases the container had two or more labels)
- 69 • Inadequate training, including training of medical gas filling personnel as well as  
70 delivery personnel
- 71 • Inadequate finished product testing
- 72 • Inadequate quality control unit
- 73 • Failure to qualify equipment prior to use (e.g., stainless steel hoses, large cryogenic  
74 containers)
- 75 • Inadequate written procedures for manufacturing, processing, testing

76

77 The Attachment, Medical Gas Mix-Ups, describes in detail some of the adverse events that the  
78 Agency has investigated, including mix-ups that have resulted in serious injury or death.

79

80 FDA can take several courses of action when a CGMP violation is found: (1) issue a warning  
81 letter; (2) seize gas-related products (including storage tanks, high-pressure cylinders, vehicles  
82 containing permanently mounted large cryogenic containers, tankers, and/or cryogenic home  
83 containers on the company's premises and trucks); (3) seek an injunction; and/or (4) initiate

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

84 prosecution. FDA may also recommend disapproval of certain government contracts with the  
85 manufacturer. FDA can also notify the Centers for Medicare & Medicaid Services (formerly the  
86 Health Care Financing Administration) of the violation. This may affect Medicare  
87 reimbursement for that company's products. FDA has issued numerous warning letters and on  
88 many occasions has successfully pursued seizure actions, injunctions, prosecutions, civil  
89 contempt actions, and inspectional warrants to enforce the CGMP regulations as they apply to  
90 medical gases.

91  
92

### **III. ORGANIZATION AND PERSONNEL**

94

#### **A. Responsibilities of the Quality Control Unit**

96

97 Medical gases are subject to the requirements in 21 CFR § 211.22 - Responsibilities of quality  
98 control unit (QCU).

99

100 Manufacturers must have a QCU with the responsibility and authority to approve or reject all  
101 product containers, closures, in-process materials, labeling, and drug products, the authority to  
102 review production records to ensure that no errors have occurred, or if errors have occurred, that  
103 they have been fully investigated. The QCU is responsible for approving or rejecting drug  
104 products manufactured, processed, packed, or held under contract by another company  
105 (§ 211.22(a)).

106

107 The QCU must have the responsibility for approving or rejecting all procedures or specifications  
108 affecting the identity, strength, quality, and purity of the drug product (§ 211.22(c)).

109

110 The responsibilities and procedures applicable to the QCU must be in writing and must be  
111 followed (§ 211.22(d)).

112

113 We recommend that the QCU perform more than a testing function, be independent of the  
114 production process, and have both quality assurance and quality control responsibilities. Ideally,  
115 the QCU would participate in and have final responsibility for all functions that could affect  
116 product quality. The corporate QCU would be responsible for reviewing and approving all  
117 written procedures, even those written by each individual location's organizational units.

118

119 We recommend that all individuals who are part of the QCU be identified in the manufacturer's  
120 operating procedures. In a well-structured and well-defined corporate structure, the QCU would  
121 be included as a separate unit. A small medical gas manufacturer can designate a single  
122 individual as the QCU.

123

124 We recommend that QCU individuals receive adequate CGMP training on a continuing basis,  
125 including quality assurance training.

126

#### **B. Personnel Qualifications**

128

129 Medical gases are subject to the requirements in § 211.25 - Personnel qualifications.



## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

130  
131 Each person engaged in the manufacturing, filling, processing, packing, or holding of a medical  
132 gas must have the education, training, and experience, or a combination thereof, to enable that  
133 person to perform the assigned functions. Training must be in the particular operations that the  
134 employee performs and in current good manufacturing practice regulations as they relate to the  
135 employee's functions. Training in the CGMP regulations must be conducted by qualified  
136 individuals on a continuing basis and with sufficient frequency to ensure that employees remain  
137 familiar with CGMP requirements applicable to them (§ 211.25(a)).  
138

139 FDA recommends that CGMP training not be conducted in one massive training session.  
140 Rather, it should be presented in smaller more manageable sessions held throughout the year, or  
141 at a minimum be held once a year. We recommend that the specific type of training received or  
142 covered, the time, and the attendance at each session be documented, and records of the training  
143 be maintained.  
144

145 Regulations at § 211.25(c) require an adequate number of qualified personnel be available to  
146 perform and supervise the manufacturing, processing, or holding of medical gases.  
147

148 Useful training information and training materials are available as shown below.  
149

- 150 • The following FDA Internet sites:

151 <http://www.fda.gov/cder/dmpq/gases.htm>

152 [www.fda.gov/cder/dmpq/cgmpregs.htm](http://www.fda.gov/cder/dmpq/cgmpregs.htm)

153 [www.fda.gov/oc/industry](http://www.fda.gov/oc/industry)

- 154 • Title 21 of the Code of Federal Regulations, Parts 210 and 211, available at:

155 [www.access.gpo.gov/nara](http://www.access.gpo.gov/nara)

- 156 • Qualified suppliers who offer CGMP training

- 157 • A qualified medical gas consultant or consulting firm

- 158 • Industry or professional associations

159  
160 The Agency recommends that each manufacturer establish and follow written training  
161 procedures for all truck drivers specific to their function, including CGMP training. Truck  
162 drivers responsible for delivery of medical gases should be trained to examine the drug label and  
163 distinguish between medical gases and industrial gases, prior to unloading a container.  
164

165 We recommend that all manufacturers who allow their drivers to connect large cryogenic  
166 containers to customer gas supply systems train their drivers in the specifics of those supply  
167 systems. We recommend cargo tanker drivers who fill medical gases into storage tanks also be  
168 trained.  
169

170 We recommend that an individual responsible for performing an odor test not have an ailment  
171 (e.g., a cold or allergies) that would adversely affect his or her sense of smell. Likewise,

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

172 employees responsible for performing the inspection for the standardized colors should be able  
173 to distinguish colors.

174

### **C. Consultants**

176

177 Medical gases are subject to the requirements in § 211.34 - Consultants.

178

179 Consultants advising on the manufacturing, processing, packing, or holding of medical gases  
180 must have sufficient education, training, and experience, or any combination thereof, to advise  
181 on the subject for which they are retained. A company must maintain records stating the name,  
182 address, and qualifications of any consultants and the type of services they provide (§ 211.34).

183

184 We recommend that consultants hired to provide assistance in achieving CGMP compliance  
185 have sufficient medical gas education, training, and/or experience.

186

187

## **IV. BUILDINGS AND FACILITIES**

189

### **A. Design and Construction**

191

192 Medical gases are subject to the requirements in § 211.42 - Design and construction features.

193

194 Any building or buildings used in the manufacture, processing, packing, or holding of a medical  
195 gas must be of a suitable size, construction, and location to facilitate cleaning, maintenance, and  
196 proper operations (§ 211.42(a)).

197

198 Buildings must have adequate space for the orderly placement of equipment and materials to  
199 prevent mix-ups and to prevent contamination (§ 211.42(b)).

200

201 Operations must be performed within specifically defined areas of adequate size. There must be  
202 separate or defined areas or other such control systems for the manufacturer's operations as are  
203 necessary to prevent contamination or mix-ups (§ 211.42(c)).

204

205 The Agency recommends that buildings be maintained in good physical condition, kept clean,  
206 and have a sufficient number of areas for organized sequential operations, such as a well-defined  
207 filling area and a well-defined quarantine area. The Agency also recommends the creation of  
208 quarantine areas to separate incoming medical gases, high-pressure cylinders, cryogenic  
209 containers, manufacturing equipment, rejected containers and closures, and the finished product.

210 No matter how large your operation, we recommend you avoid storing industrial gases and  
211 medical gases in close proximity to each other.

212

213 We also recommend that delivery vehicles have well-defined, separate areas for medical gases  
214 and industrial gases to prevent mix-ups from occurring. For example, medical and industrial  
215 gases could be separated physically in the delivery truck, or a manufacturer could use a unique  
216 identifier to distinguish medical gases from industrial gases. The Agency recommends the use of  
217 360-degree wrap-around label to identify medical gases in large cryogenic containers. If a

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

218 manufacturer applies a 360-degree wrap-around label to its large cryogenic containers, this could  
219 serve as the control system for preventing mix-ups, as long as a manufacturer has established  
220 adequate driver training, adequate written procedures, and proper stock inventory systems.

221

### **B. Security**

222

223  
224 Medical gas manufacturers are wholesale distributors who are subject to the requirements of  
225 § 205.50 - Minimum requirements for the storage and handling of prescription drugs and for the  
226 establishment and maintenance of prescription drug distribution records.

227

228 All facilities used for medical gas distribution must be secure from unauthorized entry  
229 (§ 205.50(b)(1)).

230

231 Entry into areas where medical gases are held must be limited to authorized personnel  
232 (§ 205.50(b)(1)(iii)). We recommend areas where nitrous oxide is held be especially secure.

233

234 The security requirements of § 205.50(b) apply to all facilities used for medical gas distribution.  
235 FDA interprets this regulation to include all facilities where loaded medical gas delivery trucks  
236 are parked prior to making deliveries, including at an employee's home when a loaded medical  
237 gas delivery truck is driven there and parked overnight for early morning runs.

238

239 A manufacturer could use an alarm system to secure the building and keep loading docks secure,  
240 rather than open and easily accessible.

241

242

### **V. EQUIPMENT**

243

#### **A. Equipment Cleaning and Maintenance**

244

245  
246  
247 Medical gases are subject to the requirements in § 211.67 - Equipment cleaning and  
248 maintenance.

249

250 Equipment must be cleaned, maintained, and sanitized at appropriate intervals to prevent  
251 malfunctions or contamination that would alter the safety, identity, strength, quality, or purity of  
252 the medical gas beyond the official or other established requirements (§ 211.67(a)).

253

254 Written procedures must be established and followed (§ 211.67(b), (4), (5), & (6)), including  
255 maintenance and cleaning schedules, removal or obliteration of previous batch identification,  
256 protection of clean equipment from contamination prior to use, and inspection of equipment for  
257 cleanliness immediately prior to use.

258

259 We recommend that equipment used in the manufacture of medical gas (e.g., manifolds, pigtailed,  
260 valve assemblies, hoses, and gauges) be cleaned at initial use and if exposed to a contaminant.

261 We recommend that hoses used to fill cryogenic containers have protective end caps to prevent  
262 contamination from insects, dirt, debris, and other materials. We also recommend that high-  
263 pressure cylinders exposed to the elements be provided with protective caps or some other

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

264 protective device, applied to the valve opening to prevent contamination. See related  
265 clarifications in § 211.80(b).

266  
267 We recommend that storage tanks (especially those installed at a health care facility, nursing  
268 home, or hospital), tractor trailers, rail cars, high-pressure cylinders, and cryogenic containers  
269 prior to the introduction of a medical gas be cleaned in the following circumstances: when they  
270 previously contained industrial gases; when they are first received, whether new or used; and  
271 when they are or could be, contaminated.

272

### **B. Equipment Calibration**

273

274  
275 Medical gases are subject to the requirements in § 211.68 - Automatic, mechanical, and  
276 electronic equipment.

277

278 Automatic, mechanical, or electronic equipment or other types of equipment, including  
279 computers, or related systems can be used in the manufacture, processing, packing, and holding  
280 of a drug product. If such equipment is used, it must be routinely calibrated, inspected, or  
281 checked according to a written program designed to ensure proper performance (§ 211.68(a)).  
282 Written records of those calibration checks and inspections must be maintained (§ 211.68(a)).

283

284 The Agency recommends that medical gas manufacturers use either the equipment manufacturer  
285 recommended calibration schedule or a schedule based on their own historical data. A company  
286 can reference the equipment manufacturer instruction manual in its written procedures if the  
287 manual is available for use at the manufacturing site.

288

289 We recommend that vacuum gauges undergo two calibrations. The first calibration, performed  
290 daily, would ensure that the needle on the gauge returns to zero. This check can be performed  
291 with no vacuum present, and recorded on either a batch production record or a separate log. The  
292 second calibration would ensure that vacuum gauges are calibrated based on standards  
293 established by the National Institute of Standards and Technology (NIST) on an annual basis at a  
294 minimum. Low pressure gauges and flow meters used in filling cryogenic home containers  
295 would not require calibration.

296

297 We recommend that thermometers be calibrated in accordance with manufacturer  
298 recommendations, and that the calibrations be documented in a separate log.

299

300 We also recommend that medical gas companies ensure that check valves used in a supply  
301 system to prevent the back flow of a foreign product or contaminant into the lines create a proper  
302 seal and cannot be compromised. This recommendation applies to check valves placed at  
303 various points in a supply line to protect the pump, manifold, or other equipment from over-  
304 pressurization or an undesirable back flow. Check valves do not need to be qualified if they are  
305 intended to act only as an added safety feature and do not prevent the cross contamination of  
306 gases or do not affect product identity, strength, purity, or quality.

307

### **C. Computerized Systems**

308  
309

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

310 Medical gases are subject to the requirements in § 211.68 - Automatic, mechanical, and  
311 electronic equipment.

312  
313 Appropriate controls must be exercised over computer or related systems to ensure that changes  
314 in master production and control records or other records are instituted only by authorized  
315 personnel (§ 211.68(b)). Input to and output from the computer or related system of records or  
316 data must be checked for accuracy (§ 211.68(b)). The degree and frequency of input/output  
317 verification must be based on the complexity and reliability of the computer or related system  
318 (§ 211.68(b)).

319  
320 The Agency recommends that computerized systems, including hardware and software, used in  
321 the manufacturing, processing, and holding of medical gases be validated. The depth and scope  
322 of the validation depends on the diversity, complexity, and significance of the computerized  
323 application. Commercially available software that has been qualified does not need the same  
324 level of testing as software that has been specifically developed for a company.

325  
326 The Agency recommends that computerized systems have sufficient controls to prevent  
327 unauthorized access or changes to data and to preclude omissions in data. The Agency also  
328 recommends that records be kept of any changes made to data, including who made the change,  
329 when the change was made, and the previous entry.

330  
331 We recommend that any change to computerized systems be made according to specified  
332 procedures and would be formally authorized, documented, and tested. We recommend that  
333 records of all changes, including modifications and enhancements made to hardware, software,  
334 and any other critical component of the system be kept as long as the manufacturer is still using  
335 that system.

336

337

### **VI. COMPONENTS, CONTAINERS, AND CLOSURES**

338

339

340

#### **A. General Recommendations**

341

342 Medical gases are subject to the requirements in §§ 211.80 - 211.94: Control of components and  
343 drug product containers and closures.

344

345 Manufacturers must have written procedures describing in sufficient detail the receipt,  
346 identification, storage, handling, sampling, testing, and approval or rejection of components and  
347 medical gas containers and closures (§ 211.80(a)). Containers and closures must at all times be  
348 handled and stored in a manner to prevent contamination (§ 211.80(b)).

349

350 Each medical gas container and closure, upon receipt and before acceptance, must be examined  
351 visually for appropriate labeling as to contents, container damage, and contamination  
352 (§ 211.82(a)). Containers and closures must be stored under quarantine until they have been  
353 tested or examined, as appropriate (§ 211.82(b)).

354

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

355 Medical gas containers and closures must be withheld from use until the lot has been sampled,  
356 tested, or examined, as appropriate, and released for use by the quality control unit (QCU)  
357 (§ 211.84(a)). The containers must be opened, sampled, and resealed in a manner designed to  
358 prevent contamination (§ 211.84(c)(2)). Each medical gas container and closure that is liable to  
359 contamination with filth, insect infestation, or other extraneous adulterant must be examined  
360 against established specifications for such contamination (§ 211.84(d)(5)).

361  
362 Rejected containers and closures must be identified and controlled under a quarantine system  
363 designed to prevent their use in manufacturing or processing operations for which they are  
364 unsuitable (§ 211.89).

365  
366 Medical gas containers and closures must not be reactive, additive, or absorptive so as to alter  
367 the safety, identity, strength, quality, or purity of the drug beyond the official or established  
368 requirements (§ 211.94(a)). Container closure systems must provide adequate protection against  
369 foreseeable external factors in storage and use that can cause deterioration or contamination of  
370 the drug product (§ 211.94(b)). Containers and closures must be clean (§ 211.94(c)).

371 Medical gas containers and closures are used repeatedly and therefore play a critical role in  
372 ensuring that the drug product provided to the patient has the appropriate identity, strength,  
373 quality, and purity. Containers and closures used for medical gases are integral parts of the drug  
374 delivery system. We recommend they undergo strict inspections and examinations prior to the  
375 introduction of the drug product. In addition, we advise medical gas manufacturers to determine  
376 valve assembly compatibility prior to installation on a high-pressure cylinder and during the  
377 lifetime of the valve.

378  
379 To avoid the possibility of contamination, we recommend that all high-pressure cylinders and  
380 cryogenic containers used for medical gases be dedicated to medical use only.

### *1. Prefill Inspections for Cylinders*

381  
382  
383 We recommend that the following prefill inspections be performed on each medical gas cylinder  
384 prior to the start of the filling operation. Cylinders failing any of these procedures would be  
385 quarantined to prevent their use in any subsequent filling operation. We recommend that  
386 medical gas manufacturers document all prefill inspections on a batch production record.

387  
388  
389 **Hydrostatic testing date inspection:** Hydrostatic tests offer assurance of the integrity of a  
390 cylinder. Ultrasonic inspection of steel high-pressure cylinders can be performed instead of  
391 internal visual and hydrostatic testing. We recommend that manufacturers consult U.S.  
392 Department of Transportation (DOT) requirements pertaining to hydrostatic testing of certain  
393 cylinders as appropriate (see, e.g., 49 CFR 180.209).

394  
395 **External examination:** We recommend that each cylinder be examined externally for dents, arc  
396 burns, dings, oil, grease, and other signs of damage, including fire or thermal damage, that can  
397 cause a cylinder to be unacceptable or unsafe for use. Any cylinder found to have any of these  
398 conditions would be removed from service and placed in an appropriate quarantine area until  
399 their suitability has been determined by the QCU.

400

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

401 **Venting or blowing down:** If any gas is present in a cylinder, venting or blowing down a  
402 cylinder can be performed until atmospheric pressure occurs. We recommend that cylinders  
403 containing liquid be inverted and drained.

404  
405 **Odor test:** The odor test is a very important prefill test for detecting the presence of any foreign  
406 gas or odor. Do not perform this test on carbon dioxide, nitrous oxide, toxic, or corrosive gases.  
407 If a cylinder is empty (contains no pressure), a medical gas can be introduced into the cylinder at  
408 a predetermined pressure, and an odor test can be performed on the resulting gas. Use only  
409 medical gases, as an industrial gas could contain industrial contaminants.

410  
411 ***Do not confuse this odor test with the finished product odor test conducted under § 211.165(a)***  
412 ***and required by the USP.***

413  
414 **Hammer or dead ring test:** One way to determine if a cylinder has internal corrosion is by  
415 performing a hammer or dead ring test. This test consists of lightly tapping the cylinder sidewall  
416 with a hammer-like instrument. A cylinder in good condition will make a clear bell-like ring,  
417 while a dull ring indicates possible internal corrosion. All cylinders that produce a dull ring  
418 would be quarantined until their suitability has been determined. This procedure cannot be  
419 performed on aluminum or fiber wrapped cylinders because the test would not indicate internal  
420 corrosion. A hammer test works best on empty unpressurized cylinders with a 10-year test date  
421 (stamped into the cylinder shoulder area). It is not necessary to test cylinders with a 5-year test  
422 date.

423  
424 **Valve assembly examination:** The Agency recommends that the valve assembly be appropriate  
425 for the medical gas being dispensed and be examined for debris, oil, or grease. The inspection  
426 would examine whether any of the threads on the valve or on top of the valve stem are damaged;  
427 whether the handwheel or valve stem is bent; and whether there are indications of damage,  
428 corrosion inside the valve, or excessive heat or fire damage.

429  
430 **Color code examination:** The following colors are used by the medical gas industry in the  
431 United States to aid in identifying a medical gas. We recommend manufacturers use them.

- 432
- 433 • Carbon Dioxide - gray;
  - 434 • Helium - brown;
  - 435 • Medical Air - yellow;
  - 436 • Nitrogen - black;
  - 437 • Nitrous Oxide - blue;
  - 438 • Oxygen - green; and
  - 439 • Blends of medical gases use a combination of the corresponding color for each  
440 component gas. For example, oxygen and carbon dioxide would be green and gray.

441  
442 Color coding alone cannot be relied on for identification of the medical gas; use color coding in  
443 addition to examining the product label on the cylinder.

444

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

445 **Label inspection:** We recommend that the label on the cylinder be inspected and that obsolete  
446 labels or labels containing outdated lot numbers be removed. A label on an empty cylinder does  
447 not need to be removed if it is in good condition and is identical to the label that will be used for  
448 the filled cylinder. We suggest that you ensure that cylinders bear only one manufacturer or  
449 filler's label and that you not apply new labels on top of an old label.

450  
451 **Residual gas removal:** We recommend that residual gases be removed from medical gas  
452 cylinders by means of a vacuum pump prior to filling a medical gas.

453  
454 All the above inspections can be documented on a batch production record.

### *2. Prefill Inspections for Cryogenic Home Containers*

455  
456  
457  
458 The FDA recommends that the following prefill inspections be performed on all cryogenic home  
459 containers (patient-specific containers):

- 460
- 461 • An external inspection for any signs of damage, oil, or grease that would cause the
- 462 container to be unacceptable for use
- 463 • An inspection of the inlet and outlet connection for any signs of damage, oil, or grease
- 464 • An inspection of the volume or quantity of contents gauge to ensure that it is operating
- 465 properly
- 466 • An inspection of the drug label to ensure correctness.
- 467

468 All the above inspections can be documented on a batch production record.

### *3. Dedication of Large Cryogenic Containers to Medical Use Only*

469  
470  
471  
472 To avoid the possibility of industrial contaminants, we recommend that large cryogenic  
473 containers used to contain medical gases be dedicated to medical service only.

### *4. Prefill Inspections of Large Cryogenic Containers*

474  
475  
476  
477 We recommend the following prefill inspections be performed on large cryogenic containers:

- 478
- 479 • An external examination for any signs of damage, oil, or grease that could cause the
- 480 container to be unacceptable for use
- 481 • An inspection of the inlet and outlet connections for any signs of damage, oil, or
- 482 grease and to ensure that they are the correct fittings for the corresponding medical
- 483 gas. Permanently attach all connections or fittings to the container.
- 484 • An inspection of the label for correctness.
- 485 • An examination for a 360-degree wrap-around label applied on the sidewall of the
- 486 cylinder, as close to the top portion of the container as possible, but below the top
- 487 weld seam. These labels are designed to repeat the drug product name (e.g., Medical
- 488 Oxygen) in the appropriate color around the entire container. See the "Color Code



## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

489 Examination” discussed above under section 1. *Prefill Inspections for Cylinders*; and  
490 in the Glossary under “Wrap-around” Label.

491  
492 We recommend all the above inspections be documented on a batch production record.

493  
494 5. *Prefill Inspections for Permanently Mounted Cryogenic Containers*

495  
496 We recommend that the following prefill inspections be performed on permanently mounted  
497 cryogenic containers:

- 498
- 499 • An external examination for any signs of damage, oil, or grease
  - 500 • An inspection of the inlet and outlet connections for any signs of damage, oil, or  
501 grease
  - 502 • An inspection of the product label
- 503

### **B. Retesting of Containers**

504

505  
506 Containers and closures must be retested or reexamined, as appropriate, for identity, strength,  
507 quality, and purity and approved or rejected by the QCU in accordance with § 211.84 as  
508 necessary (e.g., after storage for long periods or after exposure to air, heat or other conditions  
509 that might adversely affect the medical gas container or closure) (§ 211.87).

## **VII. PRODUCTION AND PROCESS CONTROLS**

510  
511

### **A. Written Procedures**

512  
513

514  
515  
516 Manufacturers must have written procedures for production and process controls designed to  
517 ensure that medical gases have the identity, strength, quality, and purity they purport or are  
518 represented to possess. These written procedures, including any changes, must be drafted,  
519 reviewed, and approved by the appropriate organizational units and reviewed and approved by  
520 the QCU (§ 211.100(a)).

521  
522 Written production and process control procedures must be followed in the execution of the  
523 various production and process control functions and must be documented at the time of  
524 performance (§ 211.100(b)).

525  
526 To guarantee batch uniformity and integrity of medical gases, written procedures must be  
527 established and followed that describe the in-process controls, and tests, or examinations to be  
528 conducted on appropriate samples of in-process materials of each batch. Such control  
529 procedures must be established to monitor the output and to validate the performance of those  
530 manufacturing processes that may be responsible for causing variability in the drug product (§  
531 211.110(a)).

532

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

533 The Agency recommends that the corporate QCU not allow the local QCU to establish and  
534 implement written procedures that have not been reviewed and approved by the corporate QCU.  
535

536 Written procedures provide a basis for the uniform performance of a function and a step-by-step  
537 description of how to perform a specific task, function, or operation, regardless of its size or  
538 complexity. We recommend the procedures be readily available to all employees and be read,  
539 understood, and followed by them.

540  
541 We recommend that a manufacturer or individual, especially a manufacturer filling multiple  
542 gases, have data on file demonstrating the amount of vacuum evacuation required to remove all  
543 contaminants from high-pressure cylinders. We also recommend that the manufacturer have data  
544 demonstrating that each different gas it fills would be removed by the established vacuum  
545 evacuation limit.

546  
547 We recommend that portable racks, such as those added to the main header or manifold via  
548 pigtails, be evaluated to ensure that the cylinders being filled on the portable rack are being  
549 properly vacuum evacuated and are being filled to the correct pressure, as indicated by the net  
550 content statement on the label.

551  
552 We recommend that automated filling systems (that is, systems that fill from large cryogenic  
553 containers into high pressure cylinders) be validated to provide assurance that the filling is done  
554 to the correct pressure.

### **B. Charge-in of Components**

555  
556  
557  
558 Written production and control procedures must include the following, which are designed to  
559 ensure that the medical gases produced have the identity, strength, quality, and purity they  
560 purport or are represented to possess (§ 211.101):

- 561  
562 • The batch must be formulated with the intent to provide not less than 100 percent of the  
563 labeled or established amount of active ingredient (§ 211.101(a)).  
564
- 565 • Components for medical gas manufacturing must be weighed, measured, or subdivided as  
566 appropriate (§ 211.101(b)).  
567
- 568 • Each component must be added to the batch by one person and verified by a second  
569 person (§ 211.101(d)).  
570

571 The Agency recommends that all high-pressure cylinders and cryogenic containers be filled  
572 according to the net content statement indicated on the label in accordance with section  
573 502(b)(2) of the act. This includes blends or mixtures of medical gases (i.e., multiple gases).  
574 The net content statement can be the same as the fill pressure or the service pressure. Refer to §  
575 201.51, Declaration of net quantity of contents, for further information.

#### *1. Temperature/Pressure Readings (Boyle's Law)*

576  
577  
578

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

579 A medical gas in a high-pressure cylinder increases in pressure as the temperature of the gas  
580 rises. Overfilled cylinders could reach dangerously high pressures if exposed to elevated  
581 temperatures, even if the pressure at room temperature is safe. This temperature rise can be  
582 properly compensated for during filling, so that the cylinder contents do not exceed the net  
583 content statement on the label. A temperature/pressure chart or other temperature/pressure  
584 calculation algorithms can be used to adjust the filling pressure so that the proper contents are  
585 achieved (this is usually stated as the pressure at 70°F with appropriate tolerances). We  
586 recommend that temperatures measured on the wall of a cylinder during filling operations not  
587 exceed 130°F. Before the filling is complete, the temperature and pressure reading would be  
588 recorded on the batch production record.

589  
590 To ensure that high-pressure cylinders have the correct contents as indicated on the label, the  
591 manufacturer can attach a thermometer to one cylinder per manifold-filling sequence and adjust  
592 the temperature and pressure readings according to a temperature pressure chart. We  
593 recommend that, when filling cylinders one at a time (also known as the cascade method), each  
594 cylinder have a thermometer attached to it.

595  
596 If a "+" symbol follows the hydrostatic testing date, the cylinder can be overfilled by 10 percent  
597 unless the valve is equipped with a fusible, metal-backed safety. It is critical not to overfill  
598 aluminum cylinders.

### 600 2. *Valve Assembly Leak Testing*

601  
602 The Agency recommends that a valve assembly leak test be performed during the cylinder filling  
603 operation. Each valve assembly would be tested for valve packing leaks, safety plug leaks, and  
604 other valve leaks using an appropriate leak detection solution. The test would be performed  
605 while the cylinder is under pressure with the cylinder valve open. The leak detection solution  
606 would be sprayed on and around the entire valve assembly. A leak would be indicated when  
607 bubbles appear in the solution. We recommend the solution be oxygen compatible and not  
608 contain any hydrocarbons. Solutions containing soap are not recommended because they can  
609 corrode the valve stem and can leave a residue.

610  
611 After the filling of high-pressure cylinders, and after all valves have been closed, we recommend  
612 a second valve assembly leak test be performed to detect any valve outlet leaks. If any leaks are  
613 detected, the cylinder would be removed from service and quarantined until repaired.

614  
615 The two valve assembly leak tests provide assurance that the cylinder contents do not leak out  
616 during storage or shipment, resulting in a partially filled or empty cylinder that would not  
617 contain sufficient contents for a patient.

### 618 3. *Heat of Compression*

620  
621 During the filling of high-pressure cylinders, we recommend a heat-of-compression check be  
622 performed by lightly touching the exterior of each and every cylinder. A warm cylinder  
623 indicates that the cylinder is filling properly; a cool or cold cylinder indicates that the cylinder  
624 may not be filling properly. Such a situation would be investigated.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

625  
626  
627  
628  
629  
630  
631  
632  
633  
634  
635  
636  
637  
638  
639  
640  
641  
642  
643  
644  
645  
646  
647  
648  
649  
650  
651  
652  
653  
654  
655  
656  
657  
658  
659  
660  
661  
662  
663  
664  
665  
666  
667  
668  
669

### **C. Calculation of Yield**

Actual yields and percentages of theoretical yield must be determined at the conclusion of each appropriate phase of manufacturing, processing, packing, or holding of medical gases. Such calculations must be performed by one person and independently verified by a second person (§ 211.103).

FDA recognizes that accurate inventory records and reconciliation of use are difficult to maintain for liquefied gases. Normal losses of gas occur through vaporization, the filling process, and venting and could reach 10 percent or more. The FDA does not expect the reconciliation to be 100 percent accurate. A manufacturer's procedures for reconciling the use of medical gases can include allowances for normal storage and operating losses. The procedures would include provisions for further investigation when unexplained discrepancies occur, such as losses beyond established normal levels.

## **VIII. PACKAGING AND LABELING CONTROLS**

### **A. Materials Examination and Usage**

Medical gases are subject to the requirements in § 211.122 - Materials examination and usage criteria.

There must be written procedures describing in sufficient detail the receipt, identification, storage, handling, sampling, and examination of labeling and packaging materials, and these written procedures must be followed. Labeling and packaging materials must be representatively sampled, and examined or tested upon receipt and before use in packaging or labeling of a medical gas (§ 211.122(a)).

Records must be maintained for each shipment received of each different labeling indicating receipt, examination, and whether accepted or rejected (§ 211.122(c)).

Labels for each different medical gas must be stored separately with suitable identification. Access to the storage area must be limited to authorized personnel (§ 211.122(d)).

Obsolete and outdated labels must be destroyed (§ 211.122(e)).

If cut labeling is used, labeling operations must include one of the following special control procedures (§ 211.122(g)):

- Dedication of labeling lines to each different strength of each different medical gas (§ 211.122(g)(1))
- Use of appropriate electronic or electromechanical equipment to conduct a 100 percent examination for correct labeling during or after completion of finishing operations

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 670       • Use of visual inspection to conduct a 100 percent examination for correct labeling during  
671       or after completion of finishing operations for hand-applied labeling. Such examination  
672       must be performed by one person and independently verified by a second person  
673       (§ 211.122(g)(3))

674  
675       Upon receipt from the printer, labels would be counted to verify the quantity received and would  
676       be examined to ensure correctness when compared against the master label.  
677

678       We recommend that labels be locked in a secure area with access limited to authorized  
679       personnel. Different medical gas labels would be stored separately. We recommend that  
680       industrial labels be stored in a separate area.  
681

682       It is industry practice to apply labels by hand, therefore, we recommend a second person verify  
683       the correctness of the label and document the verification. In light of recent deaths and injuries,  
684       this examination is critical to ensure that the correct label has been applied to a container of  
685       medical gas.  
686

### **B.       Labeling Control**

687  
688  
689       Medical gases are subject to the requirements in § 211.125 - Labeling issuance.  
690

691       Strict control must be exercised over labeling issued for use in medical gas labeling operations  
692       (§ 211.125(a)).  
693

694       Labeling materials issued for a batch must be carefully examined for identity and conformity to  
695       the labeling specified in the master or batch production records (§ 211.125(b)).  
696

697       Procedures must be used to reconcile the quantities of labeling issued, used, and returned, and  
698       must require evaluation of discrepancies found between the quantity of drug product finished  
699       and the quantity of labeling issued if the discrepancies are outside narrow preset limits based on  
700       historical operating data (§ 211.125(c)). However, this paragraph does not apply to the 360-  
701       degree wrap-around label that is applied to large cryogenic containers.  
702

703       The Agency recommends that all labels be issued by authorized personnel only. Before release  
704       of issued labels to an employee, we recommend a representative label be checked against the  
705       master label to ensure correctness.  
706

### **C.       Packaging and Labeling Operations**

707  
708  
709       Medical gases are subject to the requirements in § 211.130 - Packaging and labeling operations.  
710

711       There must be written procedures designed to ensure that correct labels and labeling are used for  
712       medical gases; such written procedures must be followed. These procedures must incorporate  
713       the following features (§ 211.130):  
714

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 715 • Prevention of mix-ups and cross contamination by physical or spatial separation from  
716 operations on other medical gases (§ 211.130(a))
- 717 • Identification of the medical gas with a lot or control number that permits determination  
718 of the history of the manufacture and control of the batch (§ 211.130(c))
- 719 • Examination of labeling materials for suitability and correctness before packaging  
720 operations, and documentation of such examination in the batch production record (§  
721 211.130(d))

722  
723 We recommend manufacturers consider each batch of medical gas a separate entity with unique  
724 filling procedures to help ensure that the batch is uniform and consistent. Assigning a single lot  
725 number to an entire day's production is not appropriate. Each manifold filling sequence; each  
726 uninterrupted filling sequence; and each filled cryogenic container, storage tank, and trailer  
727 would be considered a new lot and be assigned a unique lot number.

728  
729 In addition, we recommend each large cryogenic container containing liquid oxygen for delivery  
730 to patients at home, whether portable or permanently mounted in a van or a truck, be considered  
731 a lot and be assigned a unique lot number. Cryogenic home containers filled at a patient's home  
732 do not need a lot number. However, we recommend that cryogenic home containers filled on  
733 site or by a third party in advance for future delivery be given a lot number.

734  
735 For safety reasons, we recommend each medical gas container bear only one drug label  
736 containing the appropriate information. Do not place a current label on top of an obsolete label.

737  
738 In accordance with 502(b)(2) of the Act, all medical gas cylinders and cryogenic containers must  
739 bear a label with an accurate statement of the net contents. We recommend that the net contents  
740 appear on the body label or shoulder label and not on (1) a removable tag, (2) a certificate of  
741 analysis, or (3) a small separate sticker.

742  
743 If a medical gas company sells medical oxygen to emergency medical services for emergency  
744 use, the label would contain the statement:

745  
746 For emergency use only when administered by properly trained personnel for  
747 oxygen deficiency and resuscitation. For all other medical applications, Rx  
748 Only.<sup>3</sup>

749  
750 FDA would not prohibit the sale of medical oxygen with this labeling to emergency medical  
751 services (see Glossary for definition of an EMS) without a prescription.

752  
753 We recommend the labeling for large permanently mounted containers, trailers, and rail cars bear  
754 a statement consisting of "Name of the Medical Gas, Refrigerated Liquid USP or NF," such as  
755 "Oxygen Refrigerated Liquid USP."  
756

---

<sup>3</sup> See September 19, 1996, citizen petition response in docket 87P-0167.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

757 The Agency recommends the use of a 360-degree wrap-around label to identify medical gases in  
758 large cryogenic containers.

759

### **D. Drug Product Inspection**

761

762 Medical gases are subject to the requirements in § 211.134 – Drug product inspection.

763

764 Labeled products must be examined during finishing operations to provide assurance that  
765 containers in the lot have the correct label (§ 211.134(a)).

766

767 A representative sample of units must be collected at the completion of finishing operations and  
768 must be visually examined for correct labeling (§ 211.134(b)).

769

770 Results of these examinations must be recorded in the batch production or control records  
771 (§ 211.134(c)).

772

773 Only one medical gas label would appear on a cylinder or container, and the manufacturer of the  
774 medical gas would apply the label in accordance with section 502(b) of the act.

775

### **E. Expiration Dating**

776

777 Medical gases are subject to the requirements in § 211.137 - Expiration dating.

779

780 To ensure that a medical gas meets applicable standards of identity, strength, quality, and purity  
781 at the time of use, each container must bear an expiration date determined by appropriate  
782 stability testing described in § 211.166 (§ 211.137(a)).

783

784 Expiration dates must be related to any storage conditions stated on the label, as determined by  
785 stability studies described in § 211.166 (§ 211.137(b)).

786

787 Expiration dates must appear on the labeling in accordance with the requirements of § 201.17  
788 (§ 211.137(d)).

789

790 New drug products for investigational use are exempt from the requirements of this section,  
791 provided that they meet appropriate standards or specifications as demonstrated by stability  
792 studies during their use in clinical investigations (§ 211.137(g)).

793

794 The Agency recommends that high-pressure cylinders stored for long periods of time, such as  
795 those provided to patients as a backup to their oxygen concentrator, be monitored to ensure they  
796 contain the correct net contents (i.e., pressure). We recommend that companies, especially home  
797 care companies and durable medical equipment suppliers, establish and follow a written plan to  
798 periodically verify the pressure (i.e., net content) of each high-pressure cylinder stored at a  
799 patient's home and that the results be documented.

800

801

***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

802 **IX. HOLDING AND DISTRIBUTION**

803

804 **A. Warehousing Procedures**

805

806 Medical gases are subject to the requirements in § 211.142 - Warehousing procedures.

807

808 Manufacturers must develop and follow written procedures describing the warehousing of  
809 medical gases. Procedures must include (§ 211.142):

810

811 • Quarantine of medical gases before release by the QCU (§ 211.142(a))

812 • Storage of medical gases under appropriate conditions (§ 211.142(b))

813

814 The Agency recommends that separate areas be designated for the following: (1) empty  
815 containers, (2) full containers, (3) in-process containers, (4) different types of medical gases, (5)  
816 rejected containers and closures, (6) medical gases that have been released, and (7) medical  
817 gases that have not been released. We also recommend that industrial gases, containers, and  
818 equipment be stored separately from medical gases, containers, and equipment.

819

820 We recommend medical gas containers be stored under protective covering and not be subject to  
821 temperature extremes. Based on this recommendation, storage areas would be clean, dry, well  
822 ventilated, and free of combustible materials. Also all valve assemblies, hoses, and other  
823 relevant equipment would be protected from contamination such as insect infestation.

824

825 **B. Distribution Procedures and Recalls**

826

827 Medical gases are subject to the requirements in § 211.150 - Distribution procedures.

828

829 Manufacturers must establish and follow written procedures describing the distribution of  
830 medical gases (§ 211.150). They must include a system by which the distribution of each lot of  
831 the drug product can be readily determined to facilitate its recall if necessary (§ 211.150(b)).

832

833 We recommend that manufacturers have procedures to explain who would evaluate distribution  
834 information if a recall were necessary, how a recall would be initiated, who would be informed  
835 about the recall, and what would be done with the recalled product.

836

837 The Agency recommends that delivery vehicles have well-defined, separate areas for medical  
838 gases and industrial gases to prevent mix-ups from occurring. For example, medical and  
839 industrial gases can be separated physically in the delivery truck, or a manufacturer can use a  
840 unique identifier to distinguish medical gases from industrial gases. As mentioned above, the  
841 Agency recommends the use of a 360-degree wrap-around label to identify medical gases in  
842 large cryogenic containers. If a manufacturer applies a 360-degree wrap-around label to its large  
843 cryogenic containers, and the manufacturer has established adequate driver training, written  
844 procedures, and proper stock inventory systems, physical separation on a delivery vehicle is not  
845 critical.

846



***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

847 We recommend that handheld computer devices or computers used during distribution  
848 operations be validated to ensure proper performance.

849

850

851 **X. LABORATORY CONTROLS**

852

853 **A. General Controls**

854

855 Medical gases are subject to the requirements in § 211.160 - Laboratory control general  
856 requirements.

857

858 The establishment of any specifications, standards, sampling plans, test procedures, or other  
859 laboratory control mechanisms required by Subpart I of 21 CFR Part 211, including any change  
860 in such specifications, standards, sampling plans, test procedures, or other laboratory control  
861 mechanisms, must be drafted by the appropriate organizational unit and reviewed and approved  
862 by the QCU. The requirements in Subpart I of 21 CFR Part 211 must be followed and must be  
863 documented at the time of performance (§ 211.160(a)).

864

865 The Agency recommends that a manufacturer follow the specifications for the specific medical  
866 gas as described in the respective monograph of the current U.S. Pharmacopeia/National  
867 Formulary (USP/NF), or a manufacturer can establish its own specifications capable of  
868 producing equivalent or better-than-USP results.

869

870 Medical gases approved under a new drug application (NDA) or covered by an investigational  
871 new drug application (IND) would comply with the specifications established in the application.

872

873 Although a primary objective of the USP is to ensure the identity, strength, quality, and purity of  
874 a product, it is impossible to include in each monograph a test for every impurity, contaminant,  
875 or adulterant that might be present, including microbial contamination. Contaminants can arise  
876 from a change in the source of material or from a change in processing, or contaminants can be  
877 introduced from extraneous sources. We recommend that a manufacturer use tests suitable for  
878 detecting such occurrences in addition to the tests provided in the individual monograph (refer to  
879 the USP General Notices, Foreign Substances and Impurities).

880

881 In the past, deaths and injuries have resulted from adulterated products that contained  
882 contaminants or impurities that were not detected. In one example, a carbon dioxide (CO<sub>2</sub>)  
883 manufacturer in Tennessee failed to include an analysis for hydrogen cyanide in its finished  
884 product testing. As a result, the manufacturer released several large liquid batches of medical  
885 CO<sub>2</sub> that were contaminated with this deadly toxin. The source of this problem was the lack of  
886 an agreement between the supplier and the CO<sub>2</sub> manufacturer requiring notification of any  
887 change in the manufacturing process. Fortunately, the problem was discovered before any injury  
888 occurred. Our investigation found the supplier of the raw material had changed the  
889 manufacturing process, which resulted in elevated hydrogen cyanide levels. Because testing for  
890 hydrogen cyanide was not performed, an adulterated drug product was released.

891

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

892                   1.     *Sampling Plan*

893

894 We recommend that a sampling plan describe the following:

895

- 896       • How many cylinders or cryogenic containers will be tested
- 897       • When the testing will occur
- 898       • What acceptance criteria will be used for selecting samples
- 899       • What action will be taken if test results are outside established specifications

900

901                   2.     *USP Oxygen Monograph*

902

903 Medical gas manufacturers can establish their own testing specifications that meet or exceed the  
904 requirements of the USP or can use the USP specifications.

905

906 USP Testing Specifications - Specifications recommend that the potency of oxygen not be less  
907 than 99.0 percent by volume. Oxygen produced by the air liquefaction process is exempt from  
908 tests for carbon dioxide and carbon monoxide. However, if there is no documentation that the  
909 oxygen is produced by the air liquefaction process, we recommend that two additional impurity  
910 tests for carbon dioxide and for carbon monoxide be performed.

911

912 ***Note: The official method is explained below. The Agency recommends that you check***  
913 ***periodically with the USP to determine if the official method has changed or has been***  
914 ***modified.***

915

916 The ORSAT testing method uses a calibrated 100-ml buret, copper wire, and an ammonium  
917 chloride - ammonium hydroxide solution mixed together and equilibrated by agitation with the  
918 copper wire. Prior to the introduction of a sample from a pumped cylinder, a series of analyses  
919 (minimum of 3 runs) using a calibration standard would be performed to properly age the test  
920 solution and to eliminate any air bubbles that may have become trapped in the apparatus. The  
921 Agency recommends that a manufacturer not proceed with testing a filled or pumped cylinder  
922 until these analyses are completed. A 100-ml sample of the unknown gas would be drawn into  
923 the buret, agitated, and measured. An identification test, using a carbon dioxide detector tube,  
924 would be performed at the same time.

925

926 ***Note: The ammonium chloride - ammonium hydroxide solution used in this method would be***  
927 ***expected to bear an expiration date supported by appropriate stability studies.***

928

929 The USP Oxygen Monograph requires a finished drug product odor test to be performed on each  
930 container undergoing testing.

931

932 The Agency recommends that USP tests not be performed on an industrial gas in an attempt to  
933 convert it to a medical gas.

934

935 The accuracy of the USP procedure is  $\pm 0.1$  percent.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

936  
937  
938  
939  
940  
941  
942  
943  
944  
945  
946  
947  
948  
949  
950  
951  
952  
953  
954  
955  
956  
957  
958  
959  
960  
961  
962  
963  
964  
965  
966  
967  
968  
969  
970  
971  
972  
973  
974  
975  
976  
977  
978

### **3. Calibration of Instruments**

Laboratory controls must include the calibration of instruments, apparatus, gauges, and recording devices at suitable intervals in accordance with an established written program containing specific directions, schedules, limits for accuracy and precision, and provisions for remedial action in the event accuracy and/or precision limits are not met (§ 211.160(b)(4)).

Oxygen analyzers and other instruments can be calibrated at intervals specified in the instructions from the equipment manufacturer. The FDA recommends that gas manufacturers not use other medical or industrial gases as the basis for calibrating their instruments.

We recommend that standards be certified to ensure the proper level of precision and accuracy as reported on the certificate of analysis (COA).

We also recommend that each COA for a medical gas calibration standard be specific for that cylinder and provide the following information:

- Name and address of the calibration standard supplier
- Name of the product
- Lot number or unique identification number specific for each cylinder
- Analytical methodology used to assay the calibration standard
- Actual analytical results (for example, 99.9 percent nitrogen)
- The responsible person's signature and the date signed

#### **B. Testing and Release for Distribution**

Medical gases are subject to the requirements in § 211.165 - Testing and release for distribution.

For each batch of medical gas, there must be appropriate laboratory determination of satisfactory conformance to final specifications, prior to release (§ 211.165(a)).

The Agency recommends that each manufacturer determine the specific testing to be performed on any incoming medical gas and on medical gases delivered to a consignee, customer, or patient. We recommend that testing methods conform to official specifications (i.e., the USP testing methodology or a validated test procedure capable of producing equivalent or better-than-USP test method performance).

If batch results do not conform to specifications, retesting is not recommended unless a thorough investigation is performed in accordance with established written procedures.<sup>4</sup>

For high-pressure cylinders filled on a multiple outlet manifold, the Agency recommends that one or more cylinders from each manifold filling sequence be assayed for identity, strength, and

---

<sup>4</sup> A draft guidance, on *Investigating Out-of-Specification Test Results for Pharmaceutical Production* was issued on September 30, 1998.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

979 odor. For high-pressure cylinders filled individually, one cylinder per uninterrupted filling  
980 sequence can be tested for identity, strength, and odor.

### *1. Liquid-to-Liquid Filling; Oxygen Only*

984 This section pertains to the filling of liquid medical oxygen into cryogenic home containers,  
985 either at a patient's home (curbside), or on site. Due to the unique nature of this operation, the  
986 Agency recognizes that testing for conformance to final specifications prior to release is  
987 impractical. Therefore, FDA recommends using the following procedures.

#### *a. Testing of the incoming liquid oxygen*

991 In lieu of testing, the home care company (HCC) can (1) witness the testing for identity and  
992 strength of the large cryogenic container(s) performed by the supplier for each container  
993 received, (2) document that the testing has been witnessed, and (3) obtain a valid COA for each  
994 container. The employee responsible for witnessing the testing would have been trained on the  
995 analytical methodology used by their supplier. Training would be documented by the  
996 employee's company.

998 If the testing is not witnessed and the HCC chooses to rely on a valid COA, the Agency  
999 recommends that the HCC perform a specific identity test. The HCC would also periodically  
1000 verify the reliability of the supplier's analysis. This can be done by (1) visiting the supplier to  
1001 verify that the supplier is following appropriate written testing procedures, (2) observing the  
1002 supplier's analytical testing, including calibration of the analyzer, and (3) documenting that steps  
1003 1 and 2 have been taken. Alternatively, to periodically verify the reliability of the supplier's  
1004 analysis, the HCC can submit to a third party a sample from a recent delivery to be analyzed for  
1005 conformance with the USP requirements or established specifications.

1007 If an HCC does not follow the above methods or chooses to test the large cryogenic containers,  
1008 the Agency recommends that full testing on the incoming medical oxygen (each large cryogenic  
1009 container) be performed.

#### *b. Testing of an oxygen storage tank used to fill large vehicle-mounted cryogenic containers*

1014 If a new shipment of oxygen is combined in a storage tank with a previously received, tested,  
1015 and approved lot, we recommend that the manufacturer test the combined product and approve it  
1016 before use. If the storage tank is located on the company's premises and is used to fill vehicle-  
1017 mounted containers or cryogenic home containers, the Agency recommends an identity and  
1018 strength test be performed by sampling from the storage tank after each oxygen delivery and  
1019 prior to the filling of any cryogenic containers.

1021 After the storage tank has been tested, the company can forego testing a large cryogenic  
1022 container filled from the storage tank if:

- 1024 • No other storage tank is located on the premises

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 1025 • The container is dedicated to the delivery of medical oxygen for home care use only
- 1026 • The container has not been completely emptied (i.e., gaseous pressure below 15 pounds
- 1027 per square inch in gauge) and has not been out of service
- 1028 • A valid COA is received with each delivery and is maintained on file

1029  
1030  
1031

### c. Testing of cryogenic home containers

1032 It is important to exercise control over the home container during the filling operation. We  
1033 recommend that appropriate methods be developed to control situations where external  
1034 contamination may occur, such as failing to cap or cover the ends of a filling hose to prevent dirt,  
1035 debris, or insect infestation.

1036  
1037  
1038

Testing of a cryogenic home container is less of a concern if:

- 1039 • Liquid oxygen is the only liquid being filled on the premises
- 1040 • The incoming liquid oxygen is tested according to one of the methods outlined above
- 1041 under Testing of Incoming Liquid Oxygen or Testing of an Oxygen Storage Tank
- 1042 • The container is filled by the company that owns it

1043  
1044  
1045  
1046  
1047

If any other medical gas is filled on site or if the incoming liquid oxygen is not tested by one of the testing methods discussed above, we recommend all filled cryogenic home containers be tested for conformance with USP or established specifications.

1048 If cryogenic home containers are filled by another individual or another company prior to release  
1049 to the patient, we recommend that the manufacturer distributing the containers inspect each  
1050 container to ensure that a correct label including a lot number has been applied.

1051  
1052  
1053

### 2. *Liquid-To-Gas; Filling Large Cryogenic Containers*

1054 This section pertains to medical gas companies, such as welding supply companies, who fill  
1055 multiple gases, both industrial and medical. In this situation, the potential for mix-ups is  
1056 greatest. The Agency recommends that the incoming product be tested for full USP or  
1057 established specifications immediately after each lot is received. This can be done either by  
1058 taking a sample directly from the storage tank or by testing one cylinder from the first medical  
1059 filling sequence.

1060  
1061  
1062  
1063  
1064  
1065

Each filled large cryogenic container would be tested prior to release. Cryogenic containers usually contain a residual product and a commingling of new and old product would result in a new batch or lot. This new batch or lot would be analyzed and assigned a new batch or lot number. A valid COA would be provided with each cryogenic container.

1066  
1067

### 3. *Liquefied (gas on top of liquid) Compressed Gas*

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1068 The pressure in a closed container containing carbon dioxide or nitrous oxide increases as the  
1069 temperature rises. A cylinder filled at a safe pressure at normal temperatures can reach a  
1070 dangerously high pressure at high ambient temperatures. Therefore, the Agency recommends  
1071 that nitrous oxide and carbon dioxide be filled individually as liquids on a scale where the  
1072 pressure does not indicate the amount filled. Instead, we recommend these cylinders be filled  
1073 individually, and the weight not exceed 68 percent of the weight of water the cylinder will hold  
1074 at 60°F (15.6 C).

1075  
1076 The Agency recommends that one of the cylinders filled during an uninterrupted filling sequence  
1077 be tested for conformance with specifications prior to release. For both carbon dioxide and  
1078 nitrous oxide, a specific carbon dioxide identification test would be conducted concurrently with  
1079 the assay in accordance with USP monograph.

1080

### 1081 4. *Gas Mixtures*

1082

1083 If a product is a mixture of two gases, the Agency recommends that each cylinder of the blended  
1084 product be tested for the identity and strength of one of the gases, usually the active ingredient.  
1085 In addition, an identity test for the other gas would be performed on one cylinder from the  
1086 manifold filling sequence. For product mixtures containing three gases, each cylinder of the  
1087 blended product would be tested for the identity and strength of two of the gases, and one  
1088 cylinder from each manifold filling sequence would be tested for the identity of the third gas.

1089

### 1090 5. *Liquid Nitrogen*

1091

1092 An assay of the finished product using the official gas chromatographic method would not be  
1093 necessary for a manufacturer who receives shipments of medical nitrogen. However, we  
1094 recommend a manufacturer meet all of the following conditions:

1095

1096 • A valid COA is received with each delivery and the product is designated Nitrogen NF

1097

1098 • The filling system has dedicated lines, and these supply lines are traceable from the  
1099 storage tank to the filling manifold. If there is a possibility that another gas, either  
1100 industrial or medical, could be introduced and could contaminate the product, we  
1101 recommend that USP testing and a test for the absence of the contaminating gas be  
performed.

1102

1103 • Testing for the lack of oxygen (less than or equal to 1.0 percent) is performed with an  
oxygen analyzer that has been validated against the USP methodology

1104

1105 • Initially and at appropriate intervals, testing for complete specifications is recommended.  
1106 Once the reliability of the supplier is established, a manufacturer can rely upon the  
1107 supplier's certificate of analysis. Auditing the supplier's testing and manufacturing is an  
1108 additional measure that would be used to determine that the product complies with the  
1109 USP. This testing can be performed by the manufacturer, by a third party, or by a  
contract-testing laboratory.

1110

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1111 To ensure that they receive medical nitrogen, we recommend that manufacturers use suppliers  
1112 registered with FDA.

1113  
1114 In light of several reported injuries due to patient exposure to toxic compounds contained in a  
1115 supply of contaminated *industrial* nitrogen used to power surgical or dental equipment, the FDA  
1116 strongly recommends the use of *medical* nitrogen by hospitals and dentist offices, even when the  
1117 nitrogen is used for industrial purposes in those settings.<sup>5</sup>

### 1118 1119 **C. Alternate Testing Methods**

1120  
1121 The accuracy, sensitivity, specificity, and reproducibility of test methods employed by the  
1122 manufacturer must be established and documented. We recommend that such validation and  
1123 documentation be accomplished in accordance with § 211.194(a)(2).

1124  
1125 We recommend that any alternative testing method (e.g., spectrophotometer, handheld analyzers)  
1126 used to analyze a medical gas be compared against the official testing methodology.

1127  
1128 We also recommend that each medical gas manufacturer maintain a copy of the entire validation  
1129 study, including the actual data generated for each analyzer by model number that demonstrates  
1130 USP equivalence and any changes made to the analytical methodology, such as a different  
1131 column length or a different carrier gas. Validation of alternate methods can be performed in  
1132 accordance with USP Validation of Compendia Methods.

### 1133 1134 **D. Stability Testing**

1135  
1136 Medical gases are subject to the requirements in § 211.166 – Stability Testing.

1137  
1138 There must be a written testing program designed to assess the stability characteristics of  
1139 medical gases. The results of such stability testing must be used in determining appropriate  
1140 storage conditions and expiration dates. The written program must be followed and must include  
1141 (‘ 211.166(a)):

- 1142
- 1143 • Reliable, meaningful, and specific test methods (§ 211.166(a)(3))
  - 1144
  - 1145 • Testing of the medical gas in the same container-closure system as that in which the  
1146 medical gas is marketed (§ 211.166(a)(4))

1147  
1148 An adequate number of batches of each medical gas must be tested to determine an appropriate  
1149 expiration date, and a record of such data must be maintained (§ 211.166(b)).

1150  
1151 The Agency recommends that the testing program take into account the compatibility of the  
1152 valve assembly, the acceptability of the valve packing and the valve seal used, the type of  
1153 cylinder, and any other factor that can have an effect on the stability of the medical gas. Each  
1154 medical gas would be tested for stability in the exact container closure system that it is marketed

---

<sup>5</sup> Compressed Gas Association, Inc., Safety Alert (SA-6) Use of Nitrogen NF for Surgical Air Tools.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1155 in, such as steel high-pressure cylinders, aluminum high-pressure cylinders, and cryogenic  
1156 containers.

1157

### **E. Reserve Samples**

1159

1160 Reserve samples of compressed medical gases do not need to be retained (§ 211.170).

1161

1162

## **XI. RECORDS AND REPORTS**

1164

1165 Medical gases are subject to the requirements on records and reports in §§ 211.180 - 198.

1166

### **A. Record Retention**

1168

1169 Any production, control, or distribution record that is required to be maintained in compliance  
1170 with this part and is specifically associated with a batch of medical gas must be retained for at  
1171 least 1 year after the expiration date of the batch (§ 211.180(a)).

1172

1173 All records required under this part, or copies of such records, must be readily available for  
1174 authorized inspection during the retention period at the establishment where the activities  
1175 described in such records occurred (§ 211.180(c)). The records or copies thereof are subject to  
1176 photocopying or other means of reproduction as part of such an inspection (§ 211.180(c)).  
1177 Records that can be immediately retrieved from another location by computer or other electronic  
1178 means will be considered as meeting the requirements of this paragraph (§ 211.180(c)).

1179

1180 Records required under this part may be retained either as original records or as true copies such  
1181 as photocopies or other accurate reproductions of the original records (§ 211.180(d)).

1182

1183 Records can be kept on paper or electronically. Electronic records must comply with the  
1184 requirements of 21 CFR part 11.

1185

1186 Medical gas manufacturers are required to maintain a number of documents and records  
1187 including:

1188

- 1189 • Equipment cleaning and use logs (§ 211.182)
- 1190 • Computer and process validation data (§ 211.68)
- 1191 • Analyzer validation studies and data (§ 211.194)
- 1192 • Label reconciliation logs (§ 211.184)
- 1193 • Master production records (§ 211.186)
- 1194 • Batch production records (§ 211.188)
- 1195 • Analytical equipment calibration logs (§ 211.194(d))
- 1196 • Testing records (§ 211.194)
- 1197 • Stability studies (§ 211.194(e))
- 1198 • Complaint files (§ 211.198)

1199



## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1200 The Agency recommends that medical gas manufacturers also maintain training records and  
1201 certificates of analysis.

1202

### **B. Equipment Cleaning and Use Log**

1204

1205 Medical gases are subject to the requirements in § 211.182 – Equipment cleaning and use log.

1206

1207 A written record of major equipment cleaning, maintenance (except routine maintenance such as  
1208 lubrication and adjustments), and use must be included in individual equipment logs that show  
1209 the date, time, product, and lot number of each batch processed (§ 211.182). In cases where  
1210 dedicated equipment is employed, the records of cleaning, maintenance, and use must be part of  
1211 the batch record (§ 211.182). The persons performing and double-checking the cleaning and  
1212 maintenance must date and sign or initial the log indicating that the work was performed  
1213 (§ 211.182).

1214

1215 Equipment cleaning and use logs can be maintained for trailers, rail cars, and storage tanks,  
1216 especially those installed at a health care facility or a hospital.

1217

### **C. Component, Drug Product Container, Closure, and Labeling Records**

1219

1220 Medical gases are subject to the requirements in § 211.184 - Component, drug product container,  
1221 closure, and labeling records.

1222

1223 These records must include the following (§ 211.184):

1224

- 1225 • The identity and quantity of each shipment of each lot of medical labeling
- 1226 • The results of any test or examination performed (including those performed as required  
1227 by § 211.82(a), ' 211.84(d), or § 211.122(a)) and the conclusions derived therefrom
- 1228 • Documentation of the examination and review of labels and labeling for conformity with  
1229 established specifications in accordance with §§ 211.122(c) and 211.130(c)
- 1230 • The disposition of rejected medical gas containers, closures, and labeling

1231

### **D. Master Production and Control Records**

1232

1233 Medical gases are subject to the requirements in § 211.186 – Master Production and Control  
1234 Records.

1235

1236 To ensure uniformity from batch to batch, master production and control records for each  
1237 medical gas, including each batch size thereof, must be prepared, dated, and signed (full  
1238 signature, handwritten) by one person and independently checked, dated, and signed by a second  
1239 person (§ 211.186(a)). The preparation of master production and control records must be  
1240 described in a written procedure, and the written procedures must be followed (§ 211.186(a)).

1241

1242 Master production and control records must include (§ 211.186(b)):

1243

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 1244
- 1245
- 1246
- 1247
- 1248
- 1249
- 1250
- 1251
- A description of the medical gas containers and closures, and packaging materials, including a specimen or copy of each label and all other labeling signed and dated by the person or persons responsible for approval of such labeling (§ 211.186(b)(8))
  - Complete manufacturing and control instructions, sampling, and testing procedures, specifications, special notations, and precautions to be followed (§ 211.186(b)(9))

1252 The Agency recommends that a manual containing all of the above be on site and available to  
1253 personnel to ensure that individuals are able to perform their assigned functions

1254

### **E. Batch Production and Control Records**

1256

1257 Medical gases are subject to the requirements in § 211.188 – Batch Production and Control  
1258 Records.

1259

1260 Batch production and control records must be prepared for each batch of medical gas produced  
1261 and must include complete information relating to the production and control of each batch.

1262 These record must include (§ 211.188):

1263

- An accurate reproduction of the appropriate master production or control record, checked for accuracy, dated, and signed (§ 211.188(a))

1266

- Documentation that each significant step in the manufacture, processing, packing, or holding of the batch was accomplished, including (§ 211.188(b)):

1269

--Dates (§ 211.188(b)(1))

1271 --Inspection of the packaging and labeling area before and after use (§ 211.188(b)(6))

1272 --Complete labeling control records (§ 211.188(b)(8))

1273 --Description of medical gas containers and closures (§ 211.188(b)(9))

1274 --Any sampling performed (§ 211.188(b)(10))

1275 --Identification of the persons performing and directly supervising or checking each  
1276 significant step in the operations (§ 211.188(b)(11))

1277 --Any investigation made according to § 211.192 (§ 211.188(b)(12))

1278 --Results of examinations made in accordance with § 211.134 (§ 211.188(b)(13))

1279

1280 The Agency recommends that the records include documentation of the following:

1281

- Prefill inspections
- Number and size of the cylinders or cryogenic containers filled
- Filling inspections
- Post-fill inspections
- Lot number assigned
- Final temperature and pressure results
- Initials of the filler and/or analyst

1288

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- Signature of the individual who checked the entries for accuracy and completeness

Historically, the industry has used pumper's logs or filler's logs as batch production records. This is appropriate as long as the logs contain all of the relevant information. A batch production record acts as a snapshot of the actual production at the time of its performance. One batch production record would document the filling of high-pressure cylinders and a separate record would document the filling of cryogenic containers.

Based on the requirements of § 211.188, batch production records would include an item-by-item entry. A manufacturer would not use a single entry to indicate that all of the significant steps have been performed, nor a check mark or other symbol when an actual value should be recorded, such as temperature and pressure readings, purity, and identity results.

### **F. Production Record Review**

Medical gases are subject to the requirements in § 211.192 – Production record review.

All medical gas production and control records, including those for packaging and labeling, must be reviewed and approved by the QCU to determine compliance with all established, approved written procedures before a batch is released or distributed (§ 211.192). Any unexplained discrepancy or the failure of a batch to meet any of its specifications must be thoroughly investigated (§ 211.192). A written record of the investigation must be made and must include conclusions and follow-up (§ 211.192).

Any test result that is outside of the established limits would be considered an unexplained discrepancy or the failure of a batch to meet its specifications.

The Agency recommends that the release of a drug product from an air separation plant or unit (ASU) not be performed by a third-party consignee (usually known as a transporter or a trucking company). That is, the third-party consignee receiving the product would not sign as the ASU's QCU to release the product.

For ASUs where filling occurs at night, the ASU's QCU would be responsible for the release of the product, prior to distribution. For swap agreements, the manufacturer having its trailers filled would be responsible for and would have its own QCU review and approve the cleaning of any trailers that have contained industrial product, prior to filling with a medical gas.

### **G. Laboratory Records**

Medical gases are subject to the requirements in § 211.194 – Laboratory records.

Laboratory records must include complete data derived from all tests necessary to ensure compliance with established specifications and standards, including examinations and assays, as follows (§ 211.194(a)):

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 1334 • A description of the sample received for testing with identification of source (that is,  
1335 location from where the sample was obtained), quantity, lot number or other distinctive  
1336 code, date sample was taken, and date sample was received for testing (§ 211.194(a)(1))  
1337
- 1338 • A statement of each method used in the testing of the sample. The statement must  
1339 indicate the location of the data that establish that the methods used in the testing of the  
1340 sample meet proper standards of accuracy and reliability as applied to the product tested  
1341 (§ 211.194(a)(2))  
1342
- 1343 • A complete record of all data secured in the course of each test, including all graphs,  
1344 charts, and spectra from laboratory instrumentation, properly identified to show the  
1345 specific medical gas and lot tested (§ 211.194(a)(4))  
1346
- 1347 • A record of all calculations performed in connection with the test, including units of  
1348 measure, conversion factors, and equivalency factors (§ 211.194(a)(5))  
1349
- 1350 • The initials or signature of the person who performs each test and the date(s) the tests  
1351 were performed (§ 211.194(a)(7))  
1352
- 1353 • The initials or signature of a second person showing that the original records have been  
1354 reviewed for accuracy, completeness, and compliance with established standards (§  
1355 211.194(a)(8))  
1356

1357 Complete records must be maintained of any modification of an established method employed in  
1358 testing (§ 211.194(d)). Such records must include the reason for the modification and data to  
1359 verify that the modification produced results that are at least as accurate and reliable for the  
1360 material being tested as the established method (§ 211.194(b)).  
1361

1362 Complete records must be maintained of any testing and standardization of laboratory reference  
1363 standards, reagents, and standard solutions (§ 211.194(c)).  
1364

1365 Complete records must be maintained of the periodic calibration of laboratory instruments,  
1366 apparatus, gauges, and recording devices required by § 211.160(b)(4) (§ 211.194(d)).  
1367

1368 Complete records must be maintained of all stability testing performed in accordance with  
1369 § 211.166 (§ 211.194(e)).  
1370

1371 The Agency recommends that when using a handheld oxygen analyzer to perform an identity  
1372 test, the actual value obtained be recorded, and the manufacturer establish written procedures  
1373 describing an acceptable range that meets the accuracy of the analyzer.  
1374

1375 The suitability of all testing methods must be verified under actual conditions of use  
1376 (§ 211.194(a)(2)).  
1377

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1378 When testing is done by a gas chromatographic method specified in a USP monograph (such as  
1379 the assay method for Nitrogen NF), the Agency recommends the chromatographic system used  
1380 be adjusted to meet all system suitability requirements listed in the monograph. We recommend  
1381 that after tests are run to verify that the requirements have been met, the results be documented.  
1382 For monograph methods that lack specific system suitability requirements, the section on system  
1383 suitability in USP "Chromatography" can be used as a guide.

1384  
1385 When an alternative testing methodology is employed, we recommend that the methodology be  
1386 validated against an official test method and the method be carried out under substantially the  
1387 same conditions that prevailed during the validation study. If the testing environment is  
1388 substantially different, some additional on-site "spot check" tests of the method, perhaps with a  
1389 small number of standard gases, would help show that its performance has not been affected by  
1390 local conditions. For example, paramagnetic oxygen analyzers can give inaccurate readings  
1391 when used at high altitudes unless special adjustments are made. We recommend such an on-site  
1392 spot check also be made if the analyzer is installed as part of a control or alarm system. The  
1393 results of these tests would be fully documented. Certain changes made to instrumentation may  
1394 be substantive enough that they would be considered a change in the method itself; these  
1395 changes would require additional documentation of accuracy and reliability (see § 211.194(b),  
1396 above) or a new validation study.

### **H. Liquid Supply (Certificate of Analysis (COA))**

1397  
1398  
1399  
1400 The medical gas industry routinely relies on COAs to reduce the amount of finished product  
1401 testing performed. For example, if a COA lists all of the impurities tested for by a supplier, then  
1402 it would be unnecessary for a manufacturer to perform a test for the listed impurities on the  
1403 finished drug product. If no COA is received, the Agency recommends that the finished drug  
1404 product testing include all impurities listed in the USP monograph or established specifications  
1405 for each medical gas.

1406  
1407 In addition, the COA for medical oxygen usually contains the *air liquefaction statement* as  
1408 required by the USP, and as a result, it would be unnecessary for a manufacturer to test for  
1409 carbon dioxide and carbon monoxide impurities. If a manufacturer does not maintain the air  
1410 liquefaction statement for its medical oxygen, the Agency recommends that the manufacturer  
1411 perform testing for carbon dioxide and carbon monoxide impurities.

1412  
1413 We also recommend that a COA contain the following information and would accompany all  
1414 incoming deliveries of liquid medical gas:

- 1415
- 1416 • Supplier's name and complete address
  - 1417 • Name of the product (e.g., oxygen USP, carbon dioxide USP, nitrogen NF, nitrous oxide  
1418 USP, helium USP, or medical air USP)
  - 1419 • An air liquefaction statement, where appropriate
  - 1420 • Lot number or other unique identification number
  - 1421 • Actual analytical results for full USP monograph testing, (e.g., 99.5 percent oxygen)

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 1422       • Test method used to perform the analysis. If an analyzer is used, the specific model  
1423       number is indicated.  
1424       • Supplier's signature and the date  
1425       • Signature of the employee witnessing any testing at a supplier, if applicable  
1426

1427 If a company relies on a COA to reduce the amount of testing required by the USP, we  
1428 recommend the company establish the reliability of the supplier's analysis at appropriate  
1429 intervals. This can be accomplished by the manufacturer, by a third party, or by a contract-  
1430 testing laboratory.  
1431

### **I. Distribution Records**

1432       Medical gases are subject to the requirements in § 211.196 – Distribution records.  
1433

1436       Distribution records must contain the name and strength of the product and description of the  
1437       dosage form, name and address of the consignee, and date and quantity shipped (§ 211.196).  
1438

1439       A manufacturer must establish and follow written procedures that include a system whereby the  
1440       distribution of each lot of a medical gas can be determined if a recall becomes necessary, as  
1441       required in § 211.150. For compressed medical gases, distribution records are not required to  
1442       contain lot or control numbers.  
1443

### **J. Complaint Files**

1444       Medical gases are subject to the requirements in § 211.198 – Complaint files.  
1445

1448       Written procedures describing the handling of all written and oral complaints regarding a  
1449       medical gas must be established and followed (§ 211.198(a)). Such procedures must include  
1450       provisions for review by the QCU, of any complaint involving the possible failure of a medical  
1451       gas to meet any of its specifications and, for such a medical gas, a determination as to the need  
1452       for an investigation in accordance with § 211.192 (§ 211.198(a)). Such procedures must include  
1453       provisions for review to determine whether the complaint represents a serious and unexpected  
1454       adverse drug experience, which is required to be reported to the Food and Drug Administration  
1455       in accordance with § 301.305 (§ 211.198(a)).  
1456

1457       A written record of each complaint must be maintained in a file designated for medical gas  
1458       complaints (§ 211.198(b)). The file regarding such medical gas complaints must be maintained  
1459       at the establishment where the medical gas involved was manufactured, processed, or packed, or  
1460       the file may be maintained at another facility if the written records in the file are readily  
1461       available for inspection at that other facility (§ 211.198(b)).  
1462

1463       The Agency recommends that complaint records include, where known:  
1464

- 1465       • Name of the drug product  
1466       • Name and address of complainant

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

- 1467 • Name (and, where appropriate, title) and telephone number of the person submitting the
- 1468 complaint
- 1469 • Nature of the complaint
- 1470 • Date the complaint is received
- 1471 • Action initially taken, including dates and identity of person taking the action
- 1472 • Follow-up action taken
- 1473 • Response provided to the originator of complaint, including the date the response was sent
- 1474 • Final outcome regarding the issues raised by the complaint
- 1475

1476 Where an investigation under § 211.192 is conducted, the written record must include  
1477 the findings of the investigation and follow-up (211.198(b)(2)). The record or copy of  
1478 the record of the investigation must be maintained at the establishment where the investigation  
1479 occurred in accordance with § 211.180(c) (§ 211.198(b)(2)).  
1480

### **K. Reporting Deaths and Injuries**

1481  
1482 The following is intended to clarify current adverse event reporting requirements.  
1483

1484  
1485 In accordance with § 310.305, manufacturers of prescription medical gases must establish and  
1486 maintain records and must make reports to FDA of all serious, unexpected adverse drug  
1487 experiences, such as deaths or life-threatening adverse events, associated with the use of their  
1488 medical gases. More information can be obtained on FDA's web site, at  
1489 <http://www.fda.gov/medwatch/report/mfg.htm>.  
1490

1491 According to § 310.305(b) Definitions – an adverse drug experience is any adverse event  
1492 associated with the use of a drug in humans, whether or not considered drug related. This would  
1493 include problems with valves, such as valve seat combustion resulting in a release of chlorine  
1494 gas, contamination from cleaning solutions, and mix-ups that result in an adverse event to a  
1495 patient.  
1496

1497 We also encourage hospitals, nursing homes, and other health care facilities dispensing medical  
1498 gases to report serious adverse events and product problems associated with the use of those  
1499 gases. They can report adverse events directly to the medical gas manufacturer. Or, they can  
1500 report to **MedWatch**, the FDA's voluntary reporting program, in one of the following four ways:  
1501  
1502

- 1503 • Online at <http://www.accessdata.fda.gov/scripts/medwatch/>
- 1504 • By telephone at 1-800-FDA-1088
- 1505 • By FAX at 1-800-FDA-0178
- 1506 • By mail to:

1507 **MedWatch**  
1508 Food and Drug Administration (HF-2)  
1509 5600 Fishers Lane  
1510 Rockville, MD 20857-9787  
1511

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1512  
1513  
1514  
1515  
1516  
1517  
1518  
1519  
1520  
1521  
1522  
1523  
1524  
1525  
1526  
1527  
1528  
1529  
1530  
1531  
1532  
1533  
1534  
1535  
1536  
1537  
1538  
1539  
1540  
1541  
1542  
1543  
1544  
1545  
1546  
1547  
1548  
1549  
1550  
1551  
1552  
1553  
1554  
1555  
1556  
1557

### **XII. RETURNED AND SALVAGED DRUG PRODUCTS**

#### **A. Returned Drug Products**

Medical gases are subject to the requirements in § 211.204 – Returned drug products.

Returned medical gases must be identified as such and held (§ 211.204). If the conditions under which returned medical gases have been held, stored, or shipped before or during their return, or if the condition of the drug product, as a result of storage or shipping, casts doubt on the safety, identity, strength, quality or purity of the medical gas, the returned medical gas must be destroyed unless examination, testing, or other investigations prove the medical gas meets appropriate standards of safety, identity, strength, quality, or purity (§ 211.204).

#### **B. Drug Product Salvaging**

Medical gases are subject to the requirements in § 211.208 – Drug product salvaging.

Medical gases that have been subjected to improper storage conditions must not be salvaged and returned to the marketplace (§ 211.208).

### **XIII. AIR SEPARATION PLANTS OR UNITS (ASU)**

ASUs separate atmospheric air into the constituent gases of oxygen, nitrogen, and argon by using a purification process of cleaning, compressing, and cooling. ASUs are generally highly computerized and have very few employees in attendance during operations, which usually take place 24 hours a day, 7 days a week. ASUs are drug manufacturers and as such must comply with all relevant CGMP regulations.

The Agency recommends that an ASU that receives deliveries of a drug product into its storage tanks from outside sources perform finished product testing on the incoming supply, prior to accepting the delivery. Appropriate COAs would be maintained.

The Agency plans to develop and publish a separate guidance on the validation of the manufacturing process and computerized systems at ASUs.

### **XIV. STORAGE TANK INSTALLATIONS AT HEALTH CARE FACILITIES**

This section pertains to the installation of a storage tank that will contain a medical gas, usually oxygen, at a hospital, nursing home, or long-term health care facility. During the installation of a storage tank and associated equipment (i.e., equipment used for the delivery of medical gases — usually oxygen — to hospitals, nursing homes, clinics, and long-term health care facilities), following CGMP would be very important for manufacturers or individuals installing the storage tank. CGMP would also be important any time the system is exposed to a possible contaminant



## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1558 or impurity, such as installation of a new valve or piping. The company would determine the  
1559 stage of the installation where problems or contamination may occur and ensure compliance.  
1560 CGMP would be applicable to activities involving all equipment that is part of the medical gas  
1561 storage and delivery mechanism, including the storage tank that holds the drug product, all  
1562 related equipment such as piping and valves, and all other equipment up to the wall leading into  
1563 the facility.

1564

1565 For storage tank filling, we recommend a focus on the following aspects of CGMP:

1566

- 1567 • Establish a QCU and written procedures
- 1568 • Perform training for service technicians in their job functions and in CGMP
- 1569 • Qualify all equipment for medical use, including delivery vehicles and storage tanks
- 1570 • Audit contracted cleaning firms
- 1571 • Develop and follow detailed written procedures
- 1572 • Calibrate testing equipment
- 1573 • Test finished products prior to introduction of the drug product into the supply system
- 1574 • Use USP equivalent testing methodology
- 1575 • Log equipment cleaning and use, especially for storage tanks
- 1576 • Maintain batch production records
- 1577 • Provide COAs to the receiving facility with each delivery
- 1578 • Maintain documentation

1579

1580 If a third party is contracted to install a health care facility storage tank and associated  
1581 equipment, the supplier of the medical gas would determine whether the system has been  
1582 installed in accordance with CGMP. This determination would be made prior to introducing the  
1583 medical gas into the supply system and would be fully documented. The supply firm would  
1584 consider itself responsible for the actions of the third party installer.

1585

1586

### **XV. MEDICAL GAS MIX-UPS**

1587

1588  
1589 FDA has investigated a number of deaths and injuries resulting from medical gas mix-ups. In all  
1590 of these incidents, the injuries and deaths could have been prevented if the manufacturer had  
1591 followed the CGMP and industry standards.<sup>6</sup> Specific CGMP deviations noted repeatedly  
1592 included:

1593

- 1594 • § 211.100(a & b): Failure to establish and follow adequate written procedures
- 1595 • § 211.25(a): Failure to provide adequate CGMP training to all persons involved with the  
1596 handling and delivery of a medical gas
- 1597 • § 211.42(c): Inadequate storage areas on delivery vehicles for the storing of medical  
1598 gases and industrial gases

---

<sup>6</sup> See CGA standards.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1599  
1600 In particular, over the past several years, FDA has received reports from separate hospitals,  
1601 nursing homes, and clinics involving 7 deaths and 15 injuries to patients who were thought to be  
1602 receiving medical oxygen when in fact they were receiving a toxic industrial gas, (e.g., nitrogen).  
1603 The Agency recommends the following steps to help prevent similar deaths or injuries from  
1604 occurring:

- 1605  
1606 • Ensure that all employees involved in the manufacturing, processing, packaging, or  
1607 holding of medical gases have the education, training, and experience, or any  
1608 combination thereof, to enable them to perform their assigned functions.
- 1609 • Ensure that employees understand that they are handling a drug. Make sure they learn  
1610 how to examine the label on each container before delivering the container or connecting  
1611 the container to a supply system. Make sure employees know what to do if a label does  
1612 not match the invoice or the connections do not fit (e.g., possibly not accept the product  
1613 and/or notify their supervisor immediately).
- 1614 • Be aware that most fittings or connectors are permanently attached on all large cryogenic  
1615 containers used to deliver medical gases.
- 1616 • Never remove fittings and connectors. If an employee is unable to connect a container to  
1617 a supply system, he or she would contact the supplier immediately. This is especially  
1618 true for oxygen.
- 1619 • Ensure that written procedures are developed and followed. Train employees regularly  
1620 on how to perform the procedures.
- 1621 • Ensure that separate storage areas for medical and industrial gases are identified and used  
1622 on each delivery vehicle.
- 1623 • Make sure that all large cryogenic containers are dedicated to medical use and are not  
1624 used for industrial gases.
- 1625 • Ensure that all cryogenic containers have clear labeling, such as a 360-degree wrap-  
1626 around label on the sidewalls. The wrap-around label would be placed as close to the top  
1627 portion of the container as possible, but below the top weld seam, and would contain and  
1628 repeat the product name (e.g., *Medical Oxygen Medical Oxygen Medical Oxygen*) and  
1629 be the appropriate color (e.g., green for oxygen).
- 1630 • Make sure only one drug label is applied to a container. Never apply a label on top of  
1631 another label.
- 1632 • Provide each consignee (e.g., hospital, nursing home, and clinic) with a copy of FDA's  
1633 Public Health Advisory, *Guidance to Hospitals, Nursing Homes, and Other Health Care*  
1634 *Facilities*.

1635  
1636  
1637 **XVI. CARBON DIOXIDE AND HELIUM MANUFACTURERS AND WHOLESALE**  
1638 **DISTRIBUTORS**  
1639

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1640 Manufacturers of medical carbon dioxide and medical helium also are subject to CGMP  
1641 requirements. The Agency recommends that manufacturers perform process and computer  
1642 systems validation and have a written agreement with the raw material manufacturer to be  
1643 notified of any changes in the manufacturing process or the quality of the raw material. We also  
1644 recommend that manufacturers perform an initial *fingerprinting* or characterization of the  
1645 incoming raw material for any contaminants or impurities that could affect the quality, strength,  
1646 purity, or identity of the finished drug product.

1647  
1648 Carbon dioxide and helium manufacturers, as well as shippers, wholesale distributors, jobbers,  
1649 and transporters that fill these medical gases into or out of rail cars, storage tanks, trailers, and  
1650 containers are required to comply with CGMP, including the following:

- 1651
- 1652 • Process validation and computer systems validation (§ 211.68)
  - 1653 • Establishment of a QCU and written procedures (§ 211.22)
  - 1654 • In-process testing (§ 211.110)
  - 1655 • Lot numbering (§ 211.80(d))
  - 1656 • Written operating procedures (§§ 211.80(a) and 211.100)
  - 1657 • Calibration of analytical equipment (§ 211.160(b)(4))
  - 1658 • Testing of the finished product via USP or equivalent testing methodology (§ 211.165)
  - 1659 • Batch production records (§ 211.188)
  - 1660 • Maintaining documentation (§ 211.180)

1661  
1662 The Agency also recommends that carbon dioxide and helium manufacturers, as well as  
1663 shippers, wholesale distributor jobbers, and transporters that fill these medical gases into or out  
1664 of rail cars, storage tanks, trailers, and containers, do the following:

- 1665
- 1666 • Conduct training, including for CGMP
  - 1667 • Test residual gas in tankers, trailers, and rail cars prior to filling

1668  
1669 The Agency recommends that all tankers or trailers used for the delivery of carbon dioxide be  
1670 dedicated to medical use only.

1671

## 1672 **XVII. EMERGENCY MEDICAL SERVICE (EMS)**

1673  
1674 An EMS can follow this guidance to comply with CGMP when filling small high-pressure  
1675 cylinders. Given the limited nature of the operation, an EMS would emphasize:

- 1676
- 1677
  - 1678 • CGMP training
  - 1679 • Operating procedures
  - 1680 • Procedures for accurate labeling
  - 1681 • Receiving oxygen from reliable sources
  - 1682 • Performing pre-fill inspections
  - 1683 • Traceability, so that a recall can be performed if necessary
- 1684

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1685

### **XVIII. GAS-TO-GAS ADAPTERS**

1687

1688 For safety reasons, ***avoid the use of*** gas-to-gas adapters of any kind to circumvent the specific  
1689 medical gas valves and connections associated with a specific medical gas. The Agency  
1690 recommends that companies only use adapters that reduce or expand the connection size for a  
1691 specific medical gas while still maintaining the proper connection system. This practice would  
1692 be described in written procedures.

1693

1694 Adapters can be used when filling mixtures or blends. However, documented written procedures  
1695 detailing a system of checks will help prevent mix-ups or contamination. We recommend that  
1696 adapters be under strict control and be kept under limited access.

1697

### **XVIX. ALTERNATIVE APPROACHES**

1699

1700 As noted, this guidance represents FDA's current thinking on CGMP for medical gases. It does  
1701 not create or confer any rights for or on any person and does not operate to bind FDA or the  
1702 public. An alternative approach can be used if that approach satisfies the applicable statutes and  
1703 regulations. In the event you have ideas, questions, or concerns regarding an alternative  
1704 approach, we encourage you to contact the Director of the Division of Manufacturing and  
1705 Product Quality in the Office of Compliance of the Center for Drug Evaluation and Research at  
1706 FDA.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1707  
1708  
1709  
1710  
1711  
1712  
1713  
1714  
1715  
1716  
1717  
1718  
1719  
1720  
1721  
1722  
1723  
1724  
1725  
1726  
1727  
1728  
1729  
1730  
1731  
1732  
1733  
1734  
1735  
1736  
1737  
1738  
1739  
1740  
1741  
1742  
1743  
1744  
1745  
1746  
1747  
1748  
1749  
1750  
1751  
1752

### **ATTACHMENT: MEDICAL GAS MIX-UPS**

The purpose of this section is to highlight the serious consequences of failing to follow CGMP in the production and delivery of medical gases.

On December 7, 2000, a nursing home in Bellbrook, Ohio, reported the death of two patients and the injury of eight patients following a mix-up in the nursing home's oxygen supply system. The nursing home had received a shipment from their supplier of four cryogenic containers supposedly containing medical oxygen. Included in the delivery, however, was a cryogenic container of industrial nitrogen that bore two different labels. The nursing home was running low on oxygen and sent a maintenance employee to connect a new oxygen container to the oxygen supply system. The employee selected the nitrogen container and discovered, correctly, that he was unable to connect the container to the oxygen system. The employee removed a fitting from an empty oxygen container and installed it on the nitrogen container. He then connected the deadly product to the oxygen system. Several days later, two more patients died from exposure to industrial nitrogen, bringing the death toll from this one incident to four.

On December 6, 2000, an industrial nitrogen container was connected to the oxygen supply system at a Medical Center in Springerville, Arizona. The nitrogen container was properly labeled and had the correct nitrogen fitting. The supplier removed the nitrogen fitting and replaced it with an oxygen fitting. A female who had been undergoing a hysterectomy was coming off anesthesia when a ventilator alarm sounded. The anesthesiologist immediately removed the ventilator and started her on an ambu bag. The patient demonstrated no ill effects.

On July 12, 1999, a patient in a California hospital was undergoing dialysis treatment. Since he required a continuous supply of oxygen he was connected to the wall oxygen source during the procedure. Upon completion of dialysis, the oxygen connection was removed from the wall source and reattached to a portable cylinder. The cannula was attached to the patient's existing tracheostomy and the patient was transported to the Intensive Care Unit (ICU). When the patient arrived at the ICU, he was in ventricular fibrillation, became apneic and sustained a cardiac arrest. The patient died. An investigation found that the patient had been attached to a cylinder of carbon dioxide, not oxygen. This cylinder had a green top, was labeled for CO<sub>2</sub> and had the specific CO<sub>2</sub> valve.

On April 22, 1998, a hospital in Idaho discovered that a large cryogenic container of industrial nitrogen had been connected to their oxygen system supplying the operating rooms, labor and delivery rooms, and the emergency room. When the supplier's truck driver was unable to connect the incompatible nitrogen container fitting to the oxygen supply system, he used a wrench to disconnect the nitrogen fitting and replaced it with an oxygen fitting. Two patients died as a result of this medical gas mix-up.

In October 1997, a hospital in Nebraska received a shipment of large cryogenic containers which were supposed to contain medical oxygen. The shipment included one cryogenic container of industrial argon that was labeled as argon. The hospital was running low on oxygen and sent a maintenance employee to connect a new oxygen container to the oxygen supply system.

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1753 Without examining the label, the employee selected the argon container, and, discovering he was  
1754 unable to connect the container to the oxygen supply system, he removed a fitting from an empty  
1755 oxygen container, installed it on the argon container, and connected the deadly product to the  
1756 oxygen supply system. Argon was administered to a patient undergoing minor surgery. The  
1757 patient died.  
1758

1759 On December 2, 1996, a children's home located in New York reported adverse reactions  
1760 experienced by nine patients due to the inhalation of carbon dioxide. An employee of the home,  
1761 asked to attach a large cryogenic container of medical oxygen, unknowingly selected a carbon  
1762 dioxide container from their inventory. He noted that the fitting on the carbon dioxide container  
1763 was not compatible with the connector on the oxygen supply system. He removed an oxygen  
1764 fitting from an empty container, installed it on the carbon dioxide container, and attached it to  
1765 the oxygen supply system. Two patients were injured critically, four patients experienced  
1766 varying stages of respiratory distress, and three patients recovered with no lasting side effects.  
1767

1768 In March of 1996, 11 deaths were associated with contaminated oxygen delivered to a hospital  
1769 during installation of a new storage tank. A 500-gallon cryogenic container was temporarily  
1770 connected to the hospital's oxygen supply system with a 50-foot hose. An analysis of the 50-foot  
1771 hose tested positive for the presence of trichloroethylene (TCE), a standard cleaning chemical  
1772 that is very toxic to humans.  
1773

1774 In December of 1993, a home care company (HCC) that filled liquid oxygen containers  
1775 authorized an inadequately trained employee to obtain from their supplier a container (GP-45) of  
1776 medical oxygen. The supplier's employees did not accompany the HCC employee to the loading  
1777 dock to pick up the medical oxygen. The home care company's employee who failed to examine  
1778 the label selected a container of argon instead of a container of medical oxygen. The employee  
1779 loaded the container into the van and went to three patients' homes to fill their containers. When  
1780 he attempted to fill the cryogenic containers containing oxygen, the discharge line was not  
1781 compatible with the container fittings. The employee removed a fitting from an empty oxygen  
1782 container and attached it to the container containing argon, and was able to fill the patients'  
1783 containers with argon. The next day, the employee became aware of the argon mix-up and  
1784 retrieved all three containers before the patients used the gas.  
1785

1786 In July of 1986, a large welding supply company filled four gray-colored oxygen cylinders with  
1787 carbon dioxide (CO<sub>2</sub>). The cylinders were subsequently sent to a hospital and administered to  
1788 two patients undergoing surgery. One patient's death was attributed to CO<sub>2</sub> exposure; the other  
1789 patient was seriously injured. The cylinders had the proper medical oxygen label and the correct  
1790 oxygen valve. Some hospitals paint their cylinders a certain color to designate a specific unit or  
1791 room located within the hospital. In this case, the gray-colored cylinders denoted cylinders to be  
1792 delivered to the surgery rooms only.  
1793

1794 In May of 1983, a large welding supply company delivered and connected to a hospital a large  
1795 cryogenic container thought to contain medical oxygen. The gas was administered to a  
1796 premature infant, a 46-year-old male, and a 27-year-old female in three separate areas of the  
1797 hospital. All three patients died. Analysis of the container found that it contained argon instead  
1798 of oxygen, and the container bore two labels, one label read "Liquid O" while a second label on

***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1799 the opposite side of the container read "Argon"; the fill line of the container had an argon fitting;  
1800 and the discharge line had an oxygen fitting.

## Contains Nonbinding Recommendations

Draft — Not for Implementation

### GLOSSARY

1801

1802

1803 The following terms and definitions are provided to assist the reader in using this guidance  
1804 document.

1805

1806 **Cascading:** This operation pertains to gas-to-gas filling of high-pressure cylinders only, and  
1807 consists of a supply cylinder unit (usually called a *bank*) containing a group of *H* or *K*-sized  
1808 cylinders, a receiving cylinder unit, a filling manifold, and a vacuum evacuation pump. The first  
1809 supply cylinder's valve is opened and the gas flows into the smaller cylinder(s) to be filled until  
1810 equilibrium or the correct net contents is reached. If the smaller cylinder is not full and requires  
1811 additional pressure or contents, the second supply cylinder's valve is opened and the gas is  
1812 allowed to flow into the smaller cylinder. This process is repeated to the third, and fourth, etc.,  
1813 supply cylinder until the desired pressure or contents is reached in the smaller cylinder(s).  
1814 Individual cylinders in the bank are replaced sequentially as their respective pressures or  
1815 contents are diminished to levels that are ineffective for the transfilling operation.

1816

1817 **Certificate of analysis (COA):** A single document provided with each shipment of incoming  
1818 liquid medical gas that undergoes further processing (filling, transfilling). A COA contains all of  
1819 the required information that would allow the receiving manufacturer to determine if the medical  
1820 gas is acceptable. A COA can also reduce the amount of finished product testing a manufacturer  
1821 performs by allowing the manufacturer to rely on the contaminants or impurities testing  
1822 performed by the supplier and documented on the COA. Otherwise, a manufacturer would test  
1823 each finished drug product for all contaminants and impurities required by the USP or the  
1824 manufacturer's established specifications. See above for details.

1825

1826 **Cryogenic containers:** Containers used to hold a low-temperature, low-pressure liquid product  
1827 that are similar in design to an insulated thermos bottle with a vacuum between the inner and  
1828 outer container. They may be portable or permanently mounted in a vehicle, and are commonly  
1829 known as VGLs (vertical gas liquids), GPs (gas packs), or PLCs (portable liquid containers), or  
1830 HL119s, MDX 60s, 80s, and 190s. This does not include tankers, trailers, or rail cars.

1831

1832 **Cryogenic home containers:** Containers designed to hold liquid oxygen at a patient's home  
1833 under low pressure and very low temperature.

1834

1835 **Distributor:** An individual or a manufacturer that receives liquid and/or compressed gas in  
1836 labeled high-pressure cylinders or cryogenic containers and does not manipulate or apply a label  
1837 to the product. The product is then delivered to a patient or consignee.

1838

1839 **Emergency medical services (EMSs):** EMSs include fire departments, ambulance companies,  
1840 and rescue squads that are usually government-affiliated emergency services. EMSs transfill  
1841 medical oxygen for their own use (no other gases are filled on site other than compressed  
1842 breathing air) and administer medical oxygen to patients and/or victims in emergency situations.

1843

1844 **Handheld oxygen analyzers:** Oxygen analyzers that operate on the fuel cell, electrochemical  
1845 cell, galvanic cell, or polarographic principle. When properly calibrated, these analyzers provide



## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1846 a specific oxygen identification test result only. They do not have the required USP accuracy for  
1847 determining potency. We recommend they be validated.

1848  
1849 **Home care company/home respiratory care company (HCC):** Manufacturers that sell  
1850 durable medical equipment and usually supply liquid oxygen to patients at their home. They  
1851 may also fill high-pressure cylinders by means of cascading as a back up for their oxygen  
1852 concentrators.

1853  
1854 **Oxygen for environmental use:** Oxygen that meets USP specifications and is used to support  
1855 life artificially in environments that are normally deficient. This includes, but is not limited to,  
1856 space and space simulation capsules, deep submersibles, and scuba systems. This definition  
1857 excludes oxygen used in chambers or devices. This product is not to be used for inhalation or  
1858 the medical therapeutic treatment of humans or animals.

1859  
1860 **Oxygen for industrial use:** Oxygen not intended for inhalation or therapeutic treatment of  
1861 humans or animals. Because of the many contaminants and impurities associated with industrial  
1862 oxygen, industrial oxygen is not appropriate for breathing purposes.

1863  
1864 **Oxygen for aircraft use (Aviators Breathing Oxygen (ABO)):** Oxygen in fixed or portable  
1865 oxygen containers or systems intended for commercial or private aircraft use. ABO meets USP  
1866 specifications for oxygen and has special moisture and/or other limiting characteristics. We  
1867 recommend against the use of ABO for recreational inhalation or medical therapeutic treatment  
1868 of humans or animals.

1869  
1870 **Process validation:** Documented evidence that provides a high degree of assurance that a  
1871 specific process will consistently produce a product meeting its predetermined specifications and  
1872 quality attributes (see the FDA guidance, *General Principles of Process Validation*).

1873  
1874 **Storage tank or stand tank:** A large cryogenic stationary holding tank with a capacity of  
1875 several thousand to several million gallons/liters of a liquid product.

1876  
1877 **Uninterrupted filling sequence:** A single, continuous filling sequence with no breaks or  
1878 shutdowns occurring during the filling operation. This procedure uses the same personnel,  
1879 equipment, and lot of component. It does not apply to the filling of high-pressure cylinders on a  
1880 multiple outlet manifold or rack. The filling of nitrous oxide and carbon dioxide is covered by  
1881 this definition.

1882  
1883 **United States Pharmacopeia /National Formulary (USP/NF):** A reference containing a select  
1884 list of articles in the form of monographs. Included in each monograph are the standards for  
1885 determining the identity, strength, quality, and purity of the articles. Except for medical gases  
1886 approved under a new drug application or an investigational new drug application,  
1887 manufacturers can use the specifications for single medical gases described under the individual  
1888 medical gas monograph. Medical gas mixtures are not listed in the USP.

1889  
1890 **Wrap-around label:** A 360-degree label that encircles and is applied to the top of large  
1891 cryogenic containers. We recommend the lettering on the label be at least 2¾ inches high and

## ***Contains Nonbinding Recommendations***

*Draft — Not for Implementation*

1892 contain the name of the medical gas. The Agency recommends that the medical gas name be  
1893 repeated so that the name can be visible when viewed from all angles. We also recommend one  
1894 of the following: (1) the name of the medical gas (text) in the standard color for that medical gas  
1895 with a white background or (2) the background in the standard color for that medical gas with  
1896 the name of the medical gas (text) in white (See Color Code examination).