# NHANES 1999-2000 Public Release Data File (February 2003)

Lab26PP: Urinary Priority Pesticides (Non-Persistent Pesticide Metabolites), and Urinary Organophosphate Pesticides

## **Description**

# Priority Pesticides (Non-Persistent Pesticide Metabolites) and Urinary Organophosphate Pesticides

Pesticide residues and their metabolites in human tissues and fluids can be indicative of pesticide exposure and the total body burden of these pesticides. Little information is available concerning residential or household exposures to pesticides among the general population. Sufficient data do exist, however, from surveys or other focused research efforts to suggest that household exposure to certain common pesticides can be extensive and might be of significant public health concern. Pesticides of particular concern are: chlorpyrifos, 2,4-D, diazinon, permethrin, ortho-phenyl phenol, methyl parathion, and organophosphate pesticides.

## **Eligible Sample**

Participants aged 6 to 59 years.

#### **Data Collection Methods**

Urine specimens are processed, stored, and shipped to the Division of Environmental Health Laboratory Sciences, National Center for Environmental Health, Centers for Disease Control and Prevention for testing.

#### **Examination Protocol**

Detailed specimen collection and processing instructions are discussed in the <u>NHANES Laboratory/Medical Technologists Procedures Manual</u> (LPM). Vials are stored under appropriate frozen(- 20 degrees Centigrade) conditions until they are shipped to National Center for Environmental Health for testing.

## **Analytic Methodology**

**Priority Pesticides (non-persistent pesticide metabolites)** were measured using two mass spectrometric methods:

Hydrolysis, extraction, derivatization, chromatography, mass spectrometry
This method measures human metabolites of the pesticides chlorpyrifos, chlorpyrifosmethyl, naphthalene, and carbaryl in urine. The procedure involves enzymatic
hydrolysis of urine, extraction, derivatization, and column cleanup, followed by analysis
of the concentrate using capillary gas chromatography combined with tandem mass
spectrometry (GC-MS/MS). This procedure uses isotope dilution with carbon-13 labeled
internal standards for all analytes. During analysis, the analytes dissociate into charged
fragments (ions) that are specific to each analyte. Two ions from each analyte and two
ions from each carbon-13 labeled internal standard are monitored and the abundances
of each ion are measured. The ratios of these ions are used as criterion for evaluating
the data. By evaluating these analytes in urine, a measurement of the body burden of
each pesticide is obtained. For example, the amount of 1-naphthol observed in a
person's urine indicates exposure to naphthalene, a common household pesticide used
in moth repellents.

Organic Liquid Extraction, Liquid Chromatography, Tandem Mass Spectrometry This method is used for determining concentrations of select pesticide metabolites in the human population. Urine samples are prepared by extracting with an organic solvent and concentration of the solute. The concentrate is analyzed using liquid chromatography combined with tandem mass spectrometric (LC\MS\MS) determination.

The metabolites (MS parent compound) are: Atrazine mercapturate (ATZ), malathion dicarboxylic acid (MDA), 2-isopropyl-4-methyl-6-hydroxypyrimidine (IMPY), 3-phenoxybenzoic acid (3PBA), and 2,4-dichlorophenoxyacetic acid (2,4-D).

**Urinary Organophosphate** Pesticides uses Azeotropic Codistillation of Urine, Derivitization, Gas Chromatography-Tandem Mass Spectrometric method. This method can be used to quantify the six-dialkylphosphate human metabolites of organophosphate pesticides in urine. The procedure involves an azeotropic distillation of urine, derivatization using chloroiodopropane, and analysis of the chloropropyl phosphate esters using capillary gas chromatography combined with tandem mass spectrometry (GC-MS/MS). This procedure uses isotope dilution with deuterium or <sup>13</sup>C labeled internal standards for all analytes. During analysis, the protonated analytes dissociate into charged fragments (ions) that are specific to each analyte. Two ions from each analyte and two ions from each labeled internal standard are monitored and the abundances of each ion are measured. The ratios of these ions are used as a criterion for evaluating the data. By evaluating these analytes in urine, a combined measurement of the body burden of common organophosphate pesticides (cumulative exposure) is obtained.

# **Analytic Notes**

Measures of Priority Pesticides (non-persistent pesticide metabolites), and Urinary Organophosphate Pesticides are assessed in participants aged 6 to 59 years. In 1999 the sampling rate was one-fourth for participants 12 to 59 years and in 2000 the

sampling rate was one-third for participants 12 to 59 years. For participants aged 6 to 11 the sampling rate was one-half for both years.

Use the special weights included in this data file when analyzing data. Read the "Special Sample Weights for this Dataset" instructions provided before beginning analysis.

#### **Detection Limits**

For the six nonspecific organophosphate metabolites, URXOP1-URXOP6, the detection limit was constant. In cases, where the result was below the limit of detection, the value for that variable is the detection limit divided by the square root of two.

In most cases, the detection limit was variable for the other pesticide metabolites in the data set. Two variables are provided for each of these analytes. The variable named URD\_\_\_LC indicates whether the results was below the limit of detection. There are two values: "0" and "1"; "1" indicates that the result was below the limit of detection. The other variable named URX\_\_\_ provides the analytic result for that analyte. In cases, where the result was below the limit of detection, the value for that variable is the detection limit divided by the square root of two.

### **Special Sample Weights for this Dataset**

Special sample weights are required to analyze these data properly. Please check the Analytic Guidelines posted for the file on the NHANES website.

#### **Special Notes for this Dataset**

The analysis of NHANES 1999-2000 laboratory data must be conducted with the key survey design and basic demographic variables. The NHANES 1999-2000 Household Questionnaire Data Files contain demographic data, health indicators, and other related information collected during household interviews. The Household Questionnaire Data Files also contain all survey design variables and sample weights required to analyze these data. The Phlebotomy Examination file includes auxiliary information on duration of fasting, the time of day of the venipuncture, and the conditions precluding venipuncture. The Household Questionnaire and Phlebotomy Exam files may be linked to the laboratory data file using the unique survey participant identifier SEQN.